

Lonza Rockland, Inc. www.lonza.com biotechserv@lonza.com Tech Service: 800-521-0390 Customer Service: 800-638-8174 Document # 18113-0808-04 Rockland, ME 04841 USA

ProSieve[®] Protein Markers Instructions for Use

Introduction

ProSieve® Protein Markers are a novel set of 10 proteins produced by recombinant expression in bacteria. They are designed for accurate sizing of proteins in both SDS polyacrylamide gel electrophoresis and Western blotting. These recombinant proteins do not contain oligosaccharides which can cause anomalous migration, heterogeneous "fuzzy" bands, and inaccurate size estimation. Further, each protein in the set has a known, exact mass of 5, 10, 15, 25, 35, 50, 75, 100, 150 and 225 kDa. Each vial contains 550 μg of total protein with 50 μg of each individual protein except for the 50 kDa protein which has twice this amount (100 μg) and serves as a reference point.

The protein standards are supplied at their working dilution in a sample loading buffer consisting of 125 mM Tris/HCl pH 6.8, 2% SDS, 10% glycerol, 200 mM 2-mercaptoethanol, and 0.007% bromphenol blue.

Quantity: 550 μ g of protein in 500 μ l, sufficient for 100 lanes of a standard protein minigel.

Recommended Procedure

Bring the ProSieve® Protein Standards to room temperature by briefly warming them in their vial in a 37°C water bath.

For Coomassie® Brilliant Blue R-250 protein staining, load

5 μ l-10 μ l per lane of a protein minigel, and 10 μ l-30 μ l per lane for larger format gels.

For SYPRO® Red Protein Gel Stain, or silver staining, load

1 μ l-3 μ l and 2 μ l-10 μ l per lane for mini and large format gels, respectively.

Ordering Information

| Catalog No. | Description | Size |
|-------------|---------------------------|--------|
| 50547 | ProSieve® Protein Markers | 500 µl |

Related Products

PAGEr® Gold Precast Gels ProSieve® 50 Gel Solution SYPRO® Red Gel Stain SYPRO® Orange Gel Stain SYPRO® Tangerine Gel Stain AccuGENE® Tris-Glycine SDS Buffer

For Research Use Only. Not for use in diagnostic procedures.

SYPRO is a trademark of Molecular Probes, Inc. Coomassie is a trademark of Imperial Chemical Industries PLC. All other trademarks herein are marks of the Lonza Group or its affilites.

Manufactured for Lonza Rockland, Inc.

© 2007 Lonza Rockland, Inc. All rights reserved.