BioResearch



## Pluripotent Stem Cells – From Research to Therapy

Webinar

16 June 2015 / Speaker: Mr. Andrew Winner 17 June 2015 / Speaker: Dr. Heiko Bueth

# LONZC

## **Forward Statement**

Certain matters discussed in this presentation may constitute forward-looking statements. These statements are based on current expectations and estimates of Lonza Group Ltd, although Lonza Group Ltd can give no assurance that these expectations and estimates will be achieved. Investors are cautioned that all forward-looking statements involve risks and uncertainty and are qualified in their entirety. The actual results may differ materially in the future from the forward-looking statements included in this presentation due to various factors. Furthermore, except as otherwise required by law, Lonza Group Ltd disclaims any intention or obligation to update the statements contained in this presentation.



## Agenda

- Background
- Considerations for iPSC researchers with clinical aspirations
- Generation of human iPSCs from peripheral blood mononuclear cells (PBMCs)
- Xeno-free, defined culturing of human PSCs
- L7<sup>™</sup> hiPSC Reprogramming and hPSC Culture System in action – Case studies

## Agenda

## Background

- Considerations for iPSC researchers with clinical aspirations
- Generation of human iPSCs from peripheral blood mononuclear cell (PBMC)
- Xeno-free, defined culturing of human PSCs
- L7<sup>TM</sup> hiPSC Reprogramming and hPSC Culture System in action Case studies

## **Pluripo**tent Stem Cells





## The Hallmarks of Pluripotent Stem Cells



## **Charac**teristics of PSCs

PSCs possess these characteristics:

- Morphology large nuclei, refractive edge, cobblestone spacing
- Marker Expression Nanog, Tra-1-60, etc
- Karyotype Analysis Normal through many passages
- Teratoma Formation Gold standard for hPSCs
- Differentiation Potential What do your cells do?

## **Human PSC Culturing**

- What are the potential issues?
  - Every day feeding
  - Feeder culture
  - Cell recovery efficient attachment after passaging
  - Differentiation bias
  - Karyotype issues
  - Viability issues
  - Clinical aspirations



## Background

- Considerations for iPSC researchers with clinical aspirations
- Generation of human iPSCs from peripheral blood mononuclear cell (PBMC)
- Xeno-free, defined culturing of human PSCs
- L7<sup>TM</sup> hiPSC Reprogramming and hPSC Culture System in action Case studies



## **The Po**tential of Pluripotent Stem Cells **for Clin**ical Applications

#### Allogeneic Therapy (HLA-matched, HLA-null)

**Autologous Therapy** 



Several economic and technical hurdles need to be overcome to make iPSC-based therapies a reality.

## **Custom** Process Development and **Manufacturing of Therapeutic Products**



## Your Product. Our Passion.

## **Clinical**-grade iPSCs – Technical Hurdles



- cGMP-compliant
- Validated

## **Clinica**l-grade iPSCs – A Roadmap



## Lonza Generation of hiPSCs under Defined, Xenofree and cGMP-compliant Conditions





## Background

- Considerations for iPSC researchers with clinical aspirations
- Generation of human iPSCs from peripheral blood mononuclear cell (PBMC)
- Xeno-free, defined culturing of human PSCs
- L7<sup>TM</sup> The Complete System in action (case studies)



## Induced Pluripotent Stem Cells (iPSCs) Basics



## Generation of iPS Cells via Transfection of PBMCs Using Nucleofector™ Technology



L7<sup>™</sup> PBMC Reprogramming Bundle

L7<sup>™</sup> hPSC Culture System



## Lonza's Portfolio Overview

#### **iPSC** Reprogramming

#### **Cell Culture Tools**



Part	Description	Size	Part	Description	Size
CC-2702	Human PBMCs	50 M Cells	FP-5007	L7™ hPSC Media BulletKit	500 mL + 5 mL
FP-5124	L7™ PBMC Priming- Recovery Kit	24 rxn	FP-5020	L7™ hPSC Matrix	1 mg
VAXP- 3024 P3 Nu (10 Nu Ve	P3 Primary Cell 4D- Nucleofector ™ Kit (100 µl Nucleocuvette™ Vessel)	24 rxn	FP-5013	L7™ hPSC Passaging Solution	100 mL
			FP-5002	L7™ hPSC Cryosolution	50 mL

## Nucleofector<sup>™</sup> Technology – The Principle

- High transfection efficiency combined with low mortality
- DNA is directed into the nucleus giving faster gene expression



**Nucleofection Process** 



## **Compo**nents of Nucleofector<sup>™</sup> Technology



### A Unique Combination:

- Nucleofector<sup>™</sup> Device
- Specific Nucleofector<sup>™</sup> Kits
- Detailed optimized protocols
- Enabling excellent transfection performance combined with high functionality!



## iPS Generation Typically Used Human Cells Types

Human Cells*	Efficiency (pmaxGFP™ Vector)	References
Dermal fibroblasts	40-95%	Yu J <i>et al.</i> (2009) Science, 324(5928):797-801 Mehta A <i>et al.</i> (2011) Cardiovasc Res 91:577-86
Keratinocytes	40-80%	
Bone marrow CD34+ cells	66%	Chou BK et al. (2011) Cell Research, 21:518-529
Cord blood CD34+ cells	82%	Chou BK <i>et al.</i> (2011) Cell Research, 21:518-529 Mack A <i>et al.</i> (2011) PlosOne 6 (11)
Bone marrow mononuclear cells		Hu K <i>et al.</i> (2011) Blood, 117(14): e109-e119 Chou BK <i>et al.</i> (2011) Cell Research, 21:518-529
Peripheral blood mononuclear cells	30-95%	Hu K <i>et al.</i> (2011) Blood, 117(14): e109-e119 Chou BK <i>et al.</i> (2011) Cell Research, 21:518-529
Neural progenitor cells	up to 90%	
Adipose-derived stem cells	65-90%	Jia et al. (2010) Nature Methods, 7:197-199



**Blue: Lonza products** 

Black: required, not included

# **PBMC** Reprogramming Protocol Part 1

L7<sup>™</sup> Priming-Recovery Basal Medium **Priming** + Supplements Day -8\* to 0 \*In some cases a priming phase of 6 or 10 might be recommended P3 4D-Nucleofector<sup>™</sup> Kit **5**x **Nucleofection** + pCE-Oct4 + pCE-Sox/KLF Day 0 + pCE-MYC/Lin28 + pCE-mp53DD + pCXB-EBNA1 L7<sup>™</sup> Priming-Recovery Basal Medium Recovery + Supplements + Enhancer A Day 0 - 2



# **PBMC** Reprogramming Protocol **Part 2**

Blue: Lonza products Black: required, not included



## **PBMC-**derived iPSC Characterization – **Attachment and Pluripotency Markers**

Morphology (P15, D3)



Alkaline Phosphatase (P15, D4)



DAPI



TRA-1-60/Oct4



Lonza LiPSC-34B

## **PBMC-**Derived iPSC Characterization – **Karyotype and Plasmid Clearance**

#### Clear of Plasmid: LiPSC-34B (P5)

 Undetectable level of plasmid confirmed by a qPCR based assay based on the level of CT value

#### Normal Karyotype: LiPSC-34B (P11)



## L7<sup>™</sup> hPSC Culture System: The Components



- L7<sup>™</sup> hPSC Media BulletKit
- L7<sup>™</sup> hPSC Matrix
- L7<sup>TM</sup> hPSC Passaging Solution
- L7<sup>TM</sup> hPSC Cryosolution

## **Every-Other-Day Feeding – Typical**, Healthy Morphology



#### **Benefits**

- Daily maintenance of hPSCs no longer necessary
- Similar look and growth characteristics as every-day feeding



## L7<sup>™</sup> hPSC Passaging Solution Provides High Re-Plating Efficiency





Dispase



L7<sup>™</sup> hPSC Passaging Solution

## The Time Required to Expand hPSCs Using L7™ hPSC Passaging Solution is Significantly Reduced



## L7<sup>™</sup> hPSC Culture System – Unaltered Directed Differentiation Capabilities



hiPSC #6.4





Neurospheres



Cardiac Troponin MYL2 DAPI





FoxA2



Rosettes



Actin MYL2 DAPI



FoxA2 DAPI



NSCs



Cardiac Troponin MYL2 DAPI



SOX17 FoxA2

#### Lineage-specific differentiation of hPSCs

- hESCs/iPSC cultured with L7<sup>™</sup> hPSC Culture System maintain potency for lineagespecific differentiation
- Cells are not locked in undifferentiated state by L7<sup>™</sup> hPSC Culture System

## Advantages of L7™ hPSC Culture System



- Every-other-day feeding
- High re-plating efficiency
- Maintenance of pluripotency marker expression
- Normal karyotype
- Unaltered directed differentiation capabilities

# Agenda

## Background

- Considerations for iPSC researchers with clinical aspirations
- Generation of human iPSCs from peripheral blood mononuclear cell (PBMC)
- Xeno-free, defined culturing of human PSCs
- L7<sup>TM</sup> hiPSC Reprogramming and hPSC Culture System in action Case studies



## L7<sup>™</sup> hPSC Culture System for Research

■ Dr. Winston Shim evaluated the L7<sup>TM</sup> hPSC Media BulletKit, Matrix and Passaging Solution on hESCs and hiPSCs



#### Read the whitepaper

#### 33 Jun-15

## hPSC Cultures on L7<sup>™</sup> Medium to Passage 4



## L7<sup>™</sup> hPSC Passaging Solution Comparisons



- Advantages of L7<sup>™</sup> hPSC Passaging Solution
- 1) Attachment ratio post-passaging (>90% in 16 hours)
- 2) High post-detachment viability (>90%)
- 3) Generates uniform-sized aggregates
- 4) Increased split ratios (1:10)

## L7<sup>™</sup> – The Complete System The Generation of iPSCs from PBMCs Using Lonza's Protocol

Dr. Yu-Chieh Wang evaluated the Lonza protocol for generating iPSCs from a particular individual for which retro- or Sendai virus-mediated attempts had previously failed



Read the whitepaper

# Generation and Characterization of PBMC418iPS1506 line

Snapshots from the reprogramming attempt





PBMC418iPS1506 on L7™ hPSC Culture System



PBMC418iPS1506 cells positively stained with biomarkers for cellular pluripotency

# Differentiation Potential of the PBMC418iPS1506 line

PBMC418iPS1506 cells were cultured in suspension as cell aggregates for 7 days

Aggregates were plated onto human fibronectin coated plates in MelDiff medium (ScienCell Research Laboratories, Carlsbad, CA) with the appropriate cytokines and growth factors for around 30 days

Differentiated cells were maintained on MelM medium (ScienCell Research Laboratories, Carlsbad, CA)



PBMC418iPS1506 cells directed differentiated into melanocytic cells.

## Summary

- iPSC-based therapies will face different manufacturing challenges depending on the phase of development requiring new technologies and processes to overcome technical and economic hurdles.
- Lonza has created a cGMP-compliant process for the generation of clinical-grade iPSCs. From this process, L7<sup>™</sup> hiPSC Reprogramming and hPSC Culture System was created for research use!
- L7<sup>TM</sup> The Complete System advantages include fully defined and Xeno-free, easy to follow PBMC reprogramming kit and protocols, every-other-day feeding, high re-plating efficiency, pluripotent expansion and maintenance, and unaltered directed differentiation capabilities

## **Support Tools and Contact Details**

#### Our Online Databases

- Cell Database: <u>http://www.lonza.com/celldatabase</u>
- Citations: <u>http://www.lonza.com/citations</u>

#### **Technical Support**

- Scientific Support Team US: <u>scientific.support@lonza.com</u> or +1 800 638 8174
- Scientific Support Team EU + ROW: <u>scientific.support.eu@lonza.com</u> or +32 87 321 611

## **Interested in Learning More?**

#### Meet with us at ISSCR (24-27 June 2015)

- Stop by Booth B04:11
- Or join our Innovation Showcase:

"Smart Technologies for Advancing Stem Cell Research and Discoveries" Thursday, 25 June 2015, 11:30 am - 12:30 pm Room A4

More on our Innovation Showcase

- More information at:
  - www.lonza.com/ISSCR
  - www.lonza.com/L7
  - www.lonza.com/stemcells



Help your cells realize their full potential. Watch Cell's Story

## Join Our Upcoming Webinar on Live Cell Imaging

#### Cell Monitoring and Recording Remotely: The CytoSMART<sup>™</sup> System

Learn how the easy-to-use CytoSMART<sup>™</sup> System for live cell imaging enables you to take time-lapse videos and images of your cell culture without needing to manually inspect your cells.

Tuesday, 22 September 2015 11 AM PDT (Los Angeles) / 2 PM EDT (New York) <u>Register here</u>

Wednesday, 23 September 2015 9 AM BST (London) / 10 AM CEST (Berlin) / 5 PM JST (Tokyo) Register here



www.lonza.com/cytosmart

## **Thank You!**