

Isolation and Cell Culture of Primary Human Chondrocytes

Technical Reference Guide

General Remarks

This protocol is suited to isolate primary adult human chondrocytes from articular cartilage. The protocol is optimized to work on cells from different donors. However, the use of tissue with severe arthritic damages may lead to reduced performance.

Recommended Culture Media

Isolation Medium

DMEM/F-12 [1:1] [1X] [Invitrogen/Gibco; Cat. No. 31330-038] supplemented with 250 ng/ml Fungizone® Antimycotic [Invitrogen, Cat. No. 15290-026] and Penicillin/Streptomycin [Penicillin: 50 U/ml; Streptomycin: 50 µg/ml].

Culture Medium

DMEM/F-12 [1:1] [1X] [Invitrogen/Gibco Cat. No. 31330-038] supplemented with 10% FCS, 50 µg/ml 2-Phospho-L-ascorbic acid trisodium salt [Fluka, Cat. No. 49752] and Penicillin Streptomycin [Penicillin: 50 U/ml; Streptomycin: 50 µg/ml].

Additional Reagents Required

Pronase Solution

Resuspend pronase [Roche, Cat. No. 1459643] in DMEM/F-12 at a final concentration of 1.0 mg/ml and sterilize by filtration [prepare 40 ml of pronase solution for up to 15 g cartilage tissue, 60 ml are required if more than 15 g cartilage tissue will be used].

Collagenase Solution

Resuspend collagenase [Serva, Cat. No. 17465] in DMEM/F-12 Medium at a final concentration of 1mg/ml and sterilize by filtration [prepare 40 ml of collagenase solution for up to 15 g cartilage tissue, 60 ml are required if more than 15 g cartilage tissue will be used].

Collagenase / Pronase Solution

Resuspend collagenase [Serva, Cat. No. 17465] and pronase [Roche, Cat. No. 1459643] at a final concentration of 1 mg/ml each in DMEM/F-12 Medium. Pass solution through a sterile filter. Use 10 ml of this solution per 10 cm culture dish.

Cell Isolation and Cultivation of Human Chondrocytes

1. Withdraw cartilage tissue under sterile conditions and transfer it into isolation medium.
2. Cut the cartilage tissue into pieces of approximately 2 x 2 mm (preferably in a petri dish) and transfer them afterwards into a sterile 250 ml glass bottle (weigh empty bottle).
3. Wash cartilage pieces twice with PBS.
4. Add 40 ml pronase solution and shake the cartilage pieces for 30 minutes at 37°C (100 – 120 rpm).
5. Incubate the cartilage with collagenase solution for 18 to 24 hours at 37°C with slow agitation (100 – 120 rpm).
6. Filtrate the cell suspension through a 70 µm filter into 50 ml falcon tubes.
7. Centrifuge the filtered cell suspension at room temperature for 10 minutes (300 xg).
8. Discard supernatant carefully and wash cell pellets twice with PBS.
9. Resuspend cells in an appropriate volume (ca. 20 ml – 50 ml) of culture medium carefully.
10. Take an aliquot of the cell suspension (10 µl) and mix it with 90 µl trypan blue to count the cells.

Note

The digestion should be as complete as possible. Incomplete digestion will decrease the quality of the cultured chondrocytes. The digest has been performed properly if the vast majority of chondrocytes does not have a visible external matrix. In addition, most cells should be adherent 12 – 24 hours post-seeding.

Cultivation of Chondrocytes

- In order to cultivate chondrocytes in high density monolayers, 1.8 x 10⁵ cells are seeded per cm². We recommend using 10 cm culture dishes.
- Cultivate cells in high density culture for at least 3 days. You can use chondrocytes cultured for up to 10 days or longer.

Preparing Cells for Nucleofection™

Note

This incubation with pronase/collagenase is necessary to detach the cells and to remove extracellular matrix. Cells will detach quite fast, but removing the extracellular matrix takes several hours. Cells surrounded by extracellular matrix may form clumps and can be identified by their roundish shape. Singularizing cells by this long incubation with pronase/collagenase improves the Nucleofection™ Performance remarkably. You may improve this procedure by pipetting the cell suspension once per hour.

1. Wash cells once with PBS.
2. Add pronase/collagenase solution (10 ml per 10 cm culture dish) and incubate the chondrocytes for 3 – 5 hours at 37°C.
3. After the collagenase/pronase treatment, chondrocytes can easily be rinsed off the substrate.
4. Wash the cells with PBS and centrifuge (10 minutes, 300 xg) the required number of cells. The chondrocytes are now ready for Nucleofection™. Resuspend them in the appropriate Nucleofector™ Solution at cell density recommended in the Optimized Protocol.

Recommended Literature

Human Chondrocyte Culture as Models of Cartilage Specific Gene Regulation Methods in Molecular Medicine 107, 69-95; Human Cell Culture Protocols; Second Edition Humana Press Inc., Totowa, NJ, USA.

Contact Information

North America

Customer Service: 800 638 8174 (toll free)
order.us@lonza.com
Scientific Support: 800 521 0390 (toll free)
scientific.support@lonza.com

Europe

Customer Service: +32 87 321 611
order.europe@lonza.com
Scientific Support: +32 87 321 611
scientific.support.eu@lonza.com

International

Contact your local Lonza distributor
Customer Service: +1 301 898 7025
scientific.support@lonza.com

International Offices

| | |
|-----------------|----------------------------|
| Australia | +61 3 9550 0883 |
| Belgium | +32 87 321 611 |
| Brazil | +55 11 2069 8800 |
| France | 0800 91 19 81 (toll free) |
| Germany | 0800 182 52 87 (toll free) |
| India | +91 40 4123 4000 |
| Ireland | 1 800 654 253 (toll free) |
| Italy | 800 789 888 (toll free) |
| Japan | +81 3 6264 0660 |
| Luxemburg | +32 87 321 611 |
| Poland | +48 781 120 300 |
| Singapore | +65 6521 4379 |
| Spain | 900 963 298 (toll free) |
| The Netherlands | 0800 022 4525 (toll free) |
| United Kingdom | 0808 234 97 88 (toll free) |

Lonza Cologne GmbH – 50829 Cologne Germany

For research use only. Not for use in diagnostic procedures. The Nucleofector™ Technology is covered by patent and/or patent pending rights owned by the Lonza Group Ltd or its affiliates. Unless otherwise noted, all trademarks herein are marks of the Lonza Group or its affiliates. Fungizone® is a registered trademark of Invitrogen. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights.