

Clonetics™ Conditionally Immortalized Human Skeletal Muscle Cell System

Cond. Immo. SkMC – Technical Information & Instructions

Contents:

Section	Description	Page
I	Introduction	1
II	General Cell Information	2
III	Quality Control	2
IV	Cell System Components	2
V	Unpacking and Storage Instructions	2
VI	Preparation of Culture Media	3
VII	Thawing of Cells / Initiation of Culture Process	3
VIII	Maintenance	3
IX	Subculturing	4
X	Cryopreservation	5
XI	Preparation of Differentiation Medium	5
XII	Differentiation	6
XIII	Ordering Information	7
XIV	Product Warranty	7

I. Introduction

Immortalized cell lines offer the possibility of an inexhaustible supply of cells that can be used as models of animal or human tissues. However, a major limitation of current technology is that the available models are almost invariably poor representations of cells found in intact tissues. One major reason for this is that the process that allows the cells to become immortal also has an impact on the characteristics of

the cell. Lonza overcomes this limitation using one and/or two patented technologies; human telomerase reverse transcriptase (hTERT) and temperature sensitive Large-T antigen (tsSV40 LTag). The temperature dependent conditional immortalization approach utilizing the large T antigen allows the immortalization to be reserved and the cells to exhibit a range of differentiated functions. hTERT provides chromosomal stability over many cell divisions while maintaining the in vivo nature of primary cells. These technologies allow the development of cell lines from a wide range of human and animal tissues and of novel models relevant to drug discovery.

Two different cell lines were developed for human skeletal muscle cells. XM13A1 cells represent a mixed population of myoblasts isolated from a non-diabetic, 31-40 year old female and were immortalized using tsSV40 LTag. These cells form large myotubes after 7-8 days of differentiation at 37°C. XM15B1 cells represent a mixed population of myoblasts isolated from a healthy, 31-40 year old male and were immortalized using combination of tsSV40 LTag and hTERT. These cells form large myotubes after 8-10 days of differentiation at 37°C.

Myotubes derived from each cell line express markers typical of differentiation such as desmin, myogenin, GLUT4, UCP 3 and Glycogen Synthase. Functional activities found in conditionally immortalized human skeletal muscle cells include insulin-stimulated glycogen synthesis and insulin-stimulated glucose uptake.

II. General Cell Information

Cat. No.	Description	Recommended Growth Media	Cryopreserved Passage Number	Immortalizing Genes	Seeding Density Upon Thaw*	Time to Subculture
XM13A1	Skeletal Muscle	SkGM™-2 BulletKit™ Medium	Passage 6	tsSV40 LTag	30,000 cells/cm ²	7-9 days
XM15B1	Skeletal Muscle	SkGM™-2 BulletKit™ Medium	Passage 15	tsSV40 LTag + hTERT	30,000 cells/cm ²	7-9 days

*Please note that alternative seeding densities may be required for subculture and/or differentiation.

III. Quality Control

Cat. No.	Description	Cells/Vial	Viability	Maximum Productive Population Doublings	Average Doubling Time	Properties
XM13A1	Skeletal Muscle	≥800,000 cells	≥80%	30 PD	60 hrs	insulin sensitive myoblasts differentiate into myotubes
XM15B1	Skeletal Muscle	≥800,000 cells	≥80%	50 PD	50 hrs	differentiate into myotubes

All cells are performance assayed and test negative for HIV-1, mycoplasma, Hepatitis-B, Hepatitis-C, bacteria, yeast and fungi. Cell viability, morphology, cell number, and proliferative capacity are measured after recovery from cryopreservation. Clonetics™ Media are formulated for optimal growth of specific types of human cells. COAs for all media products are available upon request. Please see Section XIV (Product Warranty, Page 7) for more information on Quality Control claims and guarantees.

IV. Cell Growth System Components

(Sold Separately)

- One conditionally immortalized skeletal muscle cell product – (cryopreserved)
- One Skeletal Muscle Myoblast Cell Media BulletKit™ Medium - 500 ml
Clonetics™ SkGM™-2 BulletKit™ (Lonza Catalog No. CC-3245) contains 500 ml of Skeletal Muscle Basal Medium-2 (SkBM™-2 Medium) and the following growth supplements: human Epidermal Growth Factor (hEGF), 0.5 ml; Dexamethasone, 0.5 ml; L-glutamine, 10.0 ml; Fetal Bovine Serum (FBS), 50.0 ml; Gentamicin/Amphotericin-B (GA), 0.5 ml
- One Trypsin EDTA – 100 ml (Lonza Catalog No. CC-5012)
- One Phosphate Buffered Saline Solution (PBS) – 500 ml (Lonza Catalog No. 17-516F, or similar)

NOTE: Additional components are necessary for the differentiation and cryopreservation of these cells. Please see the corresponding selections below for more information.

V. Unpacking and Storage Instructions

1. Check all containers for leakage or breakage.
2. For cryopreserved cells: Remove cryovials from the dry ice packaging and immediately place into liquid nitrogen storage. Alternatively, thaw and use the cells immediately. If no dry ice remains, please contact Customer Service.
3. BulletKit™ Medium instructions: store basal medium at 2°-8°C and SingleQuots™ Kit at ≤-20°C in a freezer that is not self-defrosting. Once thawed, SingleQuots™ Kit should be stored at 2°-8°C and added to basal medium within 72 hours. After SingleQuots™ Kit is added to basal medium, use within 1 month. Do not re-freeze.
4. Phosphate Buffered Saline can be stored at ambient temperature. The Trypsin/EDTA Solution is sterile-filtered and then stored at -20°C until shipment. Store Trypsin/EDTA at ≤-20°C in a freezer that is not self-defrosting. Trypsin/EDTA Solution has a limited shelf life or activation at 4°C. Trypsin/EDTA may thaw during transport. If,

upon arrival, Trypsin/EDTA is thawed, store at 4°C and use within 3 days or immediately aliquot and refreeze at -20°C. Trypsin/EDTA may be thawed and refrozen once.

NOTE: To keep Trypsin/EDTA fresh and active after thawing, you may aliquot it into sterile centrifuge tubes and re-freeze at -20°C.

Using media or reagents other than what is recommended will void the cell warranty. Please contact Scientific Support if you need help selecting media and/or reagents.

VI. Preparation of Culture Media

1. Decontaminate external surfaces of all vials and the medium bottle with ethanol or isopropanol.
2. To formulate Skeletal Muscle Growth Media-2 (SkGM™-2 Medium), transfer the contents of the SkGM™-2 SingleQuots™ Kit (Catalog No. CC-3244 containing human Epidermal Growth Factor [hEGF], Dexamethasone, L-glutamine, Fetal Bovine Serum [FBS], and Gentamicin/Amphotericin-B [GA]) to SkBM™-2 Basal Medium with a pipette, and rinse each vial with medium.
3. When preparing these BulletKit™ Media, it may not be possible to recover the entire volume listed for each vial. Small losses (up to 10%) should not affect the cell growth characteristics of the supplemented medium.
4. Transfer the label provided with each kit to the basal medium bottle(s) being supplemented (avoid covering the basal medium lot # and expiration date). Use it to record the date and amount of each supplement added. After SingleQuots™ Kit is added to basal medium, store at 2°-8°C and use within 1 month. Do not freeze medium.

NOTE: If there is concern that sterility was compromised during the supplementation process, the entire newly prepared growth medium may be re-filtered with a 0.2 µm filter to assure sterility. Routine re-filtration is not recommended.

VII. Thawing of Cells / Initiation of Culture Process

NOTE: For proliferation of these cells, cells must be cultured at 33°C±1°C, 5% CO₂, 90%±2% humidity. Culturing these cells at 37°C will cause permanent and irreversible differentiation.

1. The recommended seeding density when initially thawing Conditionally Immortalized SkMC from cryopreservation is 30,000 cells/cm². One ampoule of Conditionally Immortalized SkMC containing ≥800,000 cells contains enough cells to plate at least one T-25 flask.
2. To set up culture vessels, calculate the number of vessels needed based on the recommended seeding density and the surface area of the vessels being used. Add the appropriate amount of medium to the vessels (1 ml/5 cm²) and allow the vessels to equilibrate in a 33°C±1°C, 5% CO₂, 90%±2% humidity incubator for at least 30 minutes.
3. Prior to thawing cells, add 9 ml of culture medium to a 15 ml sterile centrifuge tube and allow the vessels to equilibrate to room temperature.
4. Wipe cryovial with ethanol or isopropanol before opening. In a sterile field, briefly twist the cap a quarter turn to relieve pressure, and then retighten. Quickly thaw the cryovial in a 37°C water bath being careful not to submerge the entire vial. Watch your cryovial closely; when the last sliver of ice melts, remove it. Do not submerge it completely. Thawing the cells for longer than 2 minutes results in less than optimal results.
5. Using a micropipette, gently add the thawed cell suspension to the previously prepared 15 ml sterile centrifuge tube containing 9 ml of medium equilibrated to room temperature.
6. Centrifuge at 150 x g for 5 minutes at room temperature.
7. Carefully discard the supernatant and resuspend the pellet in 1 ml of medium using a micropipette. Dispense cells into the culture vessels set up earlier. Gently rock the culture vessel to evenly distribute the cells and return to the 33°C±1°C, 5% CO₂, 90%±2% humidity incubator.
8. Change the growth medium the day after seeding.

VIII. Maintenance

1. Change the growth medium the day after seeding and every 2-3 days thereafter.
2. Warm an appropriate amount of medium to 33°C in a sterile container. Remove the medium and replace it with the warmed, fresh medium and return the flask to the incubator.

3. Avoid repeated warming and cooling of the medium. If the entire contents are not needed for a single procedure, transfer and warm only the required volume to a sterile secondary container.

IX. Subculturing

NOTE: Lonza warrants its Clonetics™ Cells only if Lonza Subculturing Reagents are used. The recommended subculturing reagents for these cells are Trypsin/EDTA (CC-5012) and Phosphate Buffered Saline (17-516F).

The following instructions are for a 25 cm² flask. Adjust all volumes accordingly for other size flasks.

1. Subculture the cells when they are 80%-90% confluent.
2. For each 25 cm² of cells to be subcultured:
 - a. Thaw 2 ml of Trypsin/EDTA and allow to come to room temperature.
 - b. Allow 5 ml of Phosphate Buffered Saline Solution (PBS) to come to room temperature.
 - c. Remove growth medium from 4°C storage and allow warming to room temperature.
 - d. Prepare new culture vessels.
3. Subculture one flask at a time. All flasks following the first flask will be subcultured following an optimization of this protocol (explained later in this procedure), based on calculated cell count, cell viability, and seeding density.
10. At this point, tap the flask against the palm of your hand to release the majority of cells from the culture surface. If only a few cells detach, you may not have let them trypsinize long enough. Wait 30 seconds and tap again. If cells still do not detach, wait and tap every 30 seconds thereafter. This entire process should take no more than 5 minutes.
11. After cells are released, quickly transfer the detached cells to a sterile 15 ml centrifuge tube containing 4 ml of growth medium.
12. Rinse the flask with a final 1-2 ml of growth medium to collect residual cells, and add this rinse to the centrifuge tube.
13. Examine the harvested flask under the microscope to make sure the harvest was successful by looking at the number of cells left behind. This should be less than 5%.
14. Centrifuge the harvested cells at 150 x g for five minutes to pellet the cells.
 - a. Aspirate most of the supernatant, except for 100-200 µl.
 - b. Flick the cryovial with your finger to loosen the pellet.
15. Dilute the cells to a final volume of 2 to 3 ml of growth medium and note the total volume of the diluted cell suspension.
16. Determine cell count and viability using a hemacytometer and Trypan Blue or cell counter. Make a note of your cell yield for later use.
17. If necessary, dilute the suspension with growth medium to achieve the desired "cells/ml" and re-count the cells.

NOTE: The following steps must be performed in a sterile field.

4. Aspirate the medium from one culture vessel.
5. Rinse the cells with 5 ml of room temperature PBS. DO NOT forget this step. The medium contains complex proteins and calcium that neutralize the trypsin.
6. Aspirate the PBS from the flask.
7. Cover the cells with 2 ml of Trypsin/EDTA solution.
8. Place the culture vessels into a 33°C humidified incubator for 3-5 minutes. Periodically examine the cell layer microscopically and check for cell detachment.
9. Allow the trypsinization to continue until approximately 90% of the cells are rounded up.

$$\text{Total \# of Viable Cells} = \frac{\text{Total cell count} \times \text{percent viability}}{100}$$

19. The number of flasks needed depends upon cell yield and seeding density. If seeding into flasks for further proliferation, the recommended density is 10,000 cells/cm². If seeding into flasks or well plates for differentiation, the recommended density is 25,000-30,000 cells/cm². Determine the total number of flasks to inoculate by using the following equation.

$$\text{Total \# of Flasks to inoculate} = \frac{\text{Total \# of viable cells}}{\text{Growth area} \times \text{Rec. Seeding Density}}$$

NOTE: If plating cells for differentiation, plates must be gelatin coated. Please see Section XII (Differentiation, Page 6) for more information.

- Use the following equation to calculate the volume of cell suspension to seed into your flasks.

$$\text{Seeding Volume} = \frac{\text{Total volume of diluted cell suspension}}{\text{\# of flasks as determined in step 18}}$$

- Prepare flasks by labeling each flask with the passage number, cell type, and date.
- Carefully transfer growth medium to new culture vessels by adding 1 ml growth medium for every 5 cm² surface area of the flask (1 ml/5 cm²).
- After mixing the diluted cells with a 5 ml pipet to ensure a uniform suspension, dispense the calculated volume into the prepared subculture flasks.
- If not using vented caps, loosen caps of flasks. Place the new culture vessels into a 33°C±1°C, 5% CO₂, 90%±2% humidity incubator.

X. Cryopreservation

NOTE: Cryopreservation may compromise cell quality and performance.

Cryopreservation Media:

Description	Base Media	DMSO	FBS
Cond. Immo. SkMC	80% SKGM™-2	10% DMSO	10% FBS

- Prepare cryopreservation media according to the chart listed above and chill to 4°C.
- Prepare freezing vials or ampoules by labeling each with the passage number, cell type and date.
- Sterile filter cryopreservation media using a 0.2 micron filter
- Harvest and centrifuge cells according to Steps 1 to 14 of Section IX (Subculturing, Page 4).
- Resuspend cells in cold cryopreservation media at 1,000,000 cells per ml.

NOTE: Work Quickly! Once exposed to the DMSO, cells become very fragile.

- Pipet aliquots (1 ml each) into freezing vials or ampoules and seal.

- Insulate aliquots with Styrofoam or propanol freezing canister.
- Store cells at -80°C overnight.
- Within 12 to 24 hours, place cells in liquid nitrogen (-200°C) for long-term storage. Cells will be compromised by storage in -80°C.

XI. Preparation of Differentiation Medium

Differentiation Medium Components

(Sold Separately)

- DMEM:F-12 - 500 ml (Lonza Catalog No. 12-719F, or similar)
- 100x Penicillin / Streptomycin Solution (10,000 units potassium penicillin and 10,000 µg streptomycin sulfate per ml) (Lonza Catalog No. 17-602E, or similar)
- Human apo-Transferrin (Sigma Catalog No. T1147, or similar)
- 3,3',5-Triiodo-L-thyronine sodium salt (Sigma Catalog no. T6397, or similar)
- D-Biotin (Sigma Catalog No. B4639, or similar)
- Calcium Pantothenate (Sigma Catalog No. C8731, or similar)
- 1M Calcium Chloride Solution (Sigma Catalog No. 21115, or similar)
- 7.5% Bovine Serum Albumin Solution (Sigma Catalog No. A8412, or similar)
- Recombinant Human Insulin Solution (Sigma Catalog No. I9278, or similar)
- 1.0N Sodium Hydroxide Solution (Sigma Catalog No. 319511, or similar)

NOTE: For ease of preparation and storage, any of the components in the differentiation media may be first concentrated as a stock solution prior to adding to the basal media so long as the appropriate final concentration is added per the below table. For preparing stock solutions, please follow the component's manufacturer's suggestions.

- Prepare 13 µg/ml 3,3',5-Triiodo-L-thyronine sodium salt stock solution by adding 1 ml 1.0N NaOH to 0.65 mg 3,3',5-Triiodo-L-thyronine sodium salt in a 50 ml polypropylene tube. Gently swirl to dissolve, then fill tube to a final volume of 50 ml with 1.0N NaOH. Stock solution should be aliquoted and stored frozen.

Component	Amount Added	Final Concentration
DMEM:F12	500 ml	N/A
100x Pen / Strep Solution	5.000 ml	100 units potassium penicillin and 100 µg streptomycin sulfate per ml
Human apo-Transferrin	5.000 mg	10 µg/ml
3,3',5-Triiodo-L-thyronine sodium salt (13 µg/ml stock)	50 µL	0.2 nM
D-Biotin	4.0315 mg	33 µM
Calcium Pantothenate	4.0460 mg	17 µM
1 M Calcium Chloride	675 µL	1.35 mM
Recombinant Human Insulin (10 mg/ml stock)*	N/A	0.3 µg/ml

* Recombinant Human Insulin must be added immediately before use at 3 µL/10 ml of complete differentiation medium

2. Prepare differentiation medium by adding each component to the 500 ml bottle of DMEM:F12 medium using the volumes in the table above. Stir as necessary to ensure all components are completely dissolved. Store reconstituted media at 4°C for up to 30 days. **Do not complete the differentiation medium with Recombinant Human Insulin until ready to start the differentiation and feed cultures (the Recombinant Human Insulin must be used fresh).**
 - a. The first batch of “complete” differentiation medium will be prepared day 0 of the differentiation process. Add fresh Recombinant Human Insulin for each subsequent media change.
 - b. Add Recombinant Human Insulin stock solution vial to the complete differentiation media at a final concentration of 0.3 µg/ml (for example add 3 µl of 10 mg/ml Recombinant Human Insulin stock solution to every 10 ml differentiation media or 6 µl of 10 mg/ml Recombinant Human Insulin stock solution to every 20 ml differentiation media).
 - c. When changing media, do not warm up the entire bottle (i.e. for a single T-75 flask; warm up 15 ml of media).

XII. Differentiation

Differentiation Components (Sold Separately)

- One conditionally immortalized skeletal muscle cell product –(cryopreserved)

- One Differentiation Medium - 500 ml (prepared as described in Section XI [Preparation of Differentiation Medium, Page 5])
- 2% Gelatin Solution (Sigma Catalog No. G1393 or similar)
- Phosphate Buffered Saline Solution (PBS) – 500 ml (Lonza Catalog No. 17-516F, or similar)

NOTE: For differentiation of these cells, cells must be cultured at 37°C±1°C, 5% CO₂, 90%±2% humidity. Culturing these cells at 37°C will cause permanent and irreversible differentiation.

1. The recommended seeding density when plating Conditionally Immortalized SkMC for differentiation is 25,000-30,000 cells/cm².
2. To set up culture vessels, calculate the number of vessels needed based on the recommended seeding density and the surface area of the vessels being used and coat with 0.1% gelatin.
 - a. Prepare a 0.1% gelatin solution by adding 9.5 ml of Phosphate Buffered Saline Solution (PBS) to 0.5 ml of 2% Gelatin Solution.
 - b. Using pipette tips, add enough of the 0.1% gelatin solution to cover the entire culture surface area (approximately 100 µL of solution per every square centimeter of culture surface area).
 - c. Incubate the plate at 37°C for two hours, then aspirate the unbound gelatin solution using a sterile pipette.
 - d. Prior to use, thoroughly wash plates with growth medium.
3. Add the appropriate amount of growth medium to the vessels (1 ml/5 cm²) and allow the vessels to equilibrate in a 33°C±1°C, 5% CO₂, 90%±2% humidity incubator for at least 30 minutes
4. Subculture cells according to Section IX (Subculturing, Page 4).
5. Culture cells under standard culturing conditions in a 33°C±1°C, 5% CO₂, 90%±2% humidity incubator in growth medium (SkGM™-2 Medium).
6. When culture has achieved 90%-100% confluence, transfer flasks to 37°C±1°C, 5% CO₂, 90%±2% humidity incubator.
7. After 24-48 hours in 37°C incubator, remove the growth medium and replace with an equal

volume of differentiation medium containing fresh recombinant human insulin.

8. Continue to culture the cells in the differentiation medium containing fresh recombinant human insulin (replacing the differentiation medium every three days, adding fresh recombinant human insulin each time) for ~5 to 10 days, or until myotubes are observed throughout the culture. The resulting differentiated cultures can be observed to contain multinucleated (more than 3 nuclei) myotubes.
9. If the myotubes are to be used in assays that require an extended period in culture, following differentiation, remove the differentiation medium and add growth medium. For achieve optimal performance, replace the growth medium every other day to maintain the cultures for ~2 to 3 weeks post differentiation. Myotube cultures are best used by 2 weeks post-differentiation.

XIII. Ordering Information

Cryopreserved Skeletal Muscle Cells:

Cat. No.	Product	Description
XM13A1	Cond. Immo. SkMC	≥800,000 cells
XM15B1	Cond. Immo. SkMC	≥800,000 cells

Skeletal Muscle Growth Media (Sold Separately):

Cat. No.	Product	Description
CC-3245	SkGM™-2 BulletKit™ Medium	500 ml SkBM™-2 Basal Medium plus CC-3244 SingleQuots™ Kit to formulate SkGM™-2 Medium (growth medium)
CC-3246	SkBM™-2 Basal Medium	Skeletal muscle myoblast basal medium-2 (500 ml)
CC-3244	SkGM™-2 SingleQuots™ Kit	Formulates 500 ml of SkBM™-2 Basal Medium to SkGM™-2 Growth Medium; contains hEGF, 0.5 ml; Dexamethasone, 0.5 ml; L-glutamine, 10.0 ml; FBS, 50.0 ml; GA, 0.5 ml.

Subculturing Reagents (Sold Separately):

Cat. No.	Product	Description
CC-5012	Trypsin/EDTA Solution	100 ml
17-516F	PBS (1X)	Phosphate Buffered Saline Solution (1X) (500 ml)

Skeletal Muscle Differentiation Components (Sold Separately):

Cat. No.	Product	Description
12-719F	DMEM:F12	DMEM:F12, 1:1 Mixture with 3.151 g/L glucose, L-glutamine, and 15 mM HEPES (500 ml)
17-516F	PBS (1X)	Phosphate Buffered Saline Solution (1X) (500 ml)

Additional components are required for differentiation please see Section XI (Preparation of Differentiation Medium, Page 5) and Section XII (Differentiation, Page 6) for a complete listing of required components.

XIV. Product Warranty

Cultures have a finite lifespan *in vitro*.

Lonza guarantees the performance of its cells in the following manner only if Clonetics™ Media and Reagents are used exclusively and the recommend protocols are followed. The performance of cells is not guaranteed if any modifications are made to the complete cell system.

1. Conditionally Immortalized SkMC will become irreversibly differentiated if cultured at 37°C. To avoid the loss of your cells and forfeiture of your warranty, proliferation of these cells must occur at 33°C.

When placing an order or for Scientific Support, please refer to the product numbers and descriptions listed above. For a complete listing of all Clonetics™ Products, refer to the Lonza website or the current Lonza catalog. To obtain a catalog, additional information or want to speak with Scientific Support, you may contact Lonza by web, e-mail, telephone, fax or mail (See page 1 for details).

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* procedures.

WARNING: CLONETICS™ AND POIETICS™ PRODUCTS CONTAIN HUMAN SOURCE MATERIAL, TREAT AS POTENTIALLY INFECTIOUS. Each donor is tested and found non-reactive by an FDA-approved method for the presence of HIV-1, hepatitis B virus and hepatitis C virus. Where donor testing is not possible, cell products are tested for the presence of viral nucleic acid from HIV, hepatitis B virus, and hepatitis C virus. Testing cannot offer complete assurance that HIV-1, hepatitis B virus, and hepatitis C virus are absent. All human-sourced products should be handled at the biological safety level 2 to minimize exposure of potentially infectious products, as recommended in the CDC-NIH manual, [Biosafety in Microbiological and Biomedical Laboratories](#), 5th ed. If you require further information, please contact your site safety officer or Scientific Support.

All trademarks herein are marks of Lonza Group or its subsidiaries.