


Medical Devices – Quick Guide

This is a step-by-step guide detailing how to test medical devices for endotoxin.

Step 3

Soak or immerse the device(s) or flush the fluid pathway with LAL Reagent Water (LRW) that has been heated to $37 \pm 1^\circ\text{C}$, keeping the LRW in contact with the surface(s) not less than 1 hour. This can be done in a depyrogenated glass or stainless steel container. Materials may be rendered endotoxin-free by heating at 250°C for 30 minutes using a validated process. Adjust the volume for rinsing or extraction according to the size and configuration of the device.



Step 1

Select at random not more than 10 devices for testing. The number of devices tested is dependent on several factors including lot size, level of control, statistical considerations and past performance.

Note: Suitability testing for each product or product family is required. Medical devices coming from cellulose-based components may contain glucans. If glucan interference is observed, the use of a glucan-blocking reagent, such as β -G-Blocker, may be needed.

Step 4

Calculate the Endotoxin Release Limit (ERL).

The endotoxin release limit for the extraction solution can be calculated using the following equation:

$$\frac{K \times N}{V} \quad \text{where:}$$

K = allowable endotoxin per device
 20.0 EU/device
 2.15 EU/device (for cerebrospinal contact)
 N = number of devices tested
 V = total volume of extraction solution

Example:
 10 devices extracted in 1,000 mL of LRW

$$\text{ERL} = \frac{20 \text{ EU/device} \times 10 \text{ devices}}{1,000 \text{ mL}} = 0.2 \text{ EU/mL}$$

Step 2

Ensure that the sample is in its final configuration and packaging, including all component parts that make up the final medical device.

Step 5

Test the extract.

The extract should be tested with a positive product control, both in duplicate. In the event that there is interference evident, it is possible to dilute the extract further up to the Maximum Valid Dilution (MVD). The MVD can be calculated using the following equation:

$$\frac{\text{ERL}}{\lambda} \quad \text{where:}$$

ERL = calculated in step 4
 λ = the value of the lowest standard used in a quantitative assay or the labeled sensitivity of the gel clot lysate used

Example:
 Extraction solution from example in step 4 tested with the kinetic chromogenic assay where the lowest standard is 0.005 EU/mL

$$\text{MVD} = \frac{0.2 \text{ EU/mL}}{0.005 \text{ EU/mL}} = 1:40$$

www.lonza.com/pharmabiotech

www.lonza.com/lal

For more technical
guidance register to our
QC Insider™ Toolbox

www.lonza.com/qcinsider

Materials, Equipment & Documents Needed

Reagents

- Endotoxin detection kit (a variety of kits is available. Please contact your local sales representative for additional information.)
- LAL Reagent Water (LRW) (# W50-640, W50-100, W50-500, W50-1000)

Accessories

- Glass dilution tubes (# N207)
- Glass reaction tubes (# N201, N205), for use with the PYROGENT™ Gel Clot Assays
- Individually wrapped serological pipettes (optional)
- Tips
- 96-well plates (# 25-340), if using a plate-based method
- Reagent reservoirs (# 00190035), if using a plate-based method
- β-G-Blocker (# N190), optional


Use pyrogen-free accessories that have been qualified for endotoxin testing.

Equipment and Software

- Pipettors
- Timer
- Vortex mixer
- Test tube rack
- For use with the PYROGENT™ Gel Clot Assay: heating block or non-circulating hot water bath
- For use with the QCL-1000™ Endpoint Chromogenic Assay: heating block with adapter (# 25-038A) and spectro- or filterphotometer with 405-410 nm filter or microplate reader
- For use with the Kinetic-QCL™ or PYROGENT™-5000 Kinetic LAL Assays: incubating absorbance microplate reader
- For use with the PyroGene™ rFC Assay: incubating fluorescence microplate reader

Supporting Documents

- Certificate of Analysis (CoA)
- Endotoxin Detection Kit Package Insert
- Endotoxin Detection Assay Procedure Quick Guide
- Endotoxin Assay Selection Guide
- United States Pharmacopeial Convention. General Chapter <161> Medical Devices – Bacterial Endotoxin and Pyrogen Tests
- United States Pharmacopeial Convention. General Chapter <85> Bacterial Endotoxins Test
- European Pharmacopoeia (EP). European Directorate for the Quality of Medicines. Chapter 2.6.14 Bacterial Endotoxins



North America

Customer Service: 800 638 8174 (toll free)
order.us@lonza.com

Scientific Support: 800 521 0390 (toll free)
scientific.support@lonza.com

Europe

Customer Service: +32 87 321 611
order.europe@lonza.com

Scientific Support: +32 87 321 611
scientific.support.eu@lonza.com

International

Contact your local Lonza distributor

Customer Service: +1 301 898 7025

Fax: +1 301 845 8291

scientific.support@lonza.com

Lonza Walkersville, Inc. – Walkersville, MD 21793

All trademarks belong to Lonza or its affiliates or to their respective third party owners. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights.

©2016 Lonza. All rights reserved.
RT-DS024 03/16
