



# The PYROSTAR™ ES-F Series



### FUJiFILM Value from Innovation

## PYROSTAR™ ES-F Series

### Unique Endotoxin-Specific LAL Reagents that combine the simplicity of the Gel-Clot method with the quantitation of a Kinetic-Turbidimetric Assav

Why Use An Endotoxin-Specific LAL Reagent?

Untreated LAL reagents react not only with endotoxin, but also with  $(1\rightarrow 3)$   $\beta$ -D-glucan, a fungal cell wall component that initiates the clotting cascade by activating the Factor G pathway. Glucan activation results in false positive endotoxin results! To eliminate this (1→3) β-D-glucan activation, various endotoxin-specific reagents have been developed by either removing Factor G or inhibiting its activation.

The reagents in this series are available in either a singletest format with reaction vials containing pre-dispensed LAL reagent or a multi-test format requiring an aliquot from the dissolved reagent to be dispensed into endotoxin-free glass test tubes. The singletest format is ideal for an assay with a few samples, and the multi-test format for a larger number of samples.





### PYROSTAR™ ES-F Series

The PYROSTAR™ ES-F series of reagents are endotoxin-specific and unreactive to (1→3) β-D-glucan. In addition, these products are formulated to be used as either a gel-clot reagent or a kinetic-turbidimetric reagent.

The reagents in this series are available in either a single-test format with reaction vials already containing pre-dispensed LAL reagent for a single measurement, or a multi-test format requiring an aliquot from the dissolved reagent to be dispensed into endotoxin-free glass test tubes. The single-test format is ideal for an assay with a few samples, and the multi-test format for a larger number of samples.

The multi-test reagent is used by dispensing 0.1 mL of dissolved LAL reagent into reaction tubes, which is then mixed after having 0.1 mL of the sample added. The single-test reagent vial is used by adding 0.2 mL of the sample directly to the reaction vial with pre-dispensed, lyophilized LAL reagent.

## PYROSTAR™ ES-F Series FUJIFILM



### 1) Gel-Clot Assay

Gel-clot LAL reagents are supplied with a specified lysate sensitivity or lambda (λ). Lysate sensitivity is the minimum concentration of endotoxin required to generate a gel-clot. Samples, standards or controls are mixed with the LAL reagent in an endotoxin-free glass test tube and incubated undisturbed at 37 ± 1 °C for 60 ± 2 minutes. Upon completion of the incubation period, the tubes are slowly inverted 180°. If a gel has formed and maintains its integrity without deformation or collapse throughout the 180° inversion, the result is scored a positive. The reaction mixture it is scored a negative if no gel has formed and it remains a liquid. Typically, a series of sample dilutions (usually 2-fold) are assayed. The highest positive dilution of each replicate is referred to as the endpoint. The endotoxin concentration in the sample is calculated by multiplying the geometric mean dilution factor by the labeled sensitivity.

The PYROSTAR™ ES-F series is available with lysate sensitivities ranging from 0.015 to 0.25 EU/mL, and are available in both single and multi-test vials.

### 2) Kinetic-Turbidimetric Assay (KTA)

The kinetic-turbidimetric assay uses the rate of change in reagent turbidity to quantitate the amount of endotoxin in samples. The onset of turbidity is referred to as the gelation time (Tg). Tg is typically measured in minutes.

The Tg is inversely related to the endotoxin concentration. The higher the endotoxin concentration, the faster the reaction mixture becomes turbid and the shorter the Tg. The lower the endotoxin concentration, the slower the reaction mixture becomes turbid the longer the Tg.

Kinetic-Turbidimetric LAL reagents are supplied with a quantitative range of endotoxin concentrations over which the amount of endotoxin in a sample can be quantitated. The Tg of standards containing a known amount of endotoxin is used to construct a standard curve relating Tg to its corresponding endotoxin concentration. The endotoxin concentration in samples and controls is calculated from their respective Tg's and the standard curve.

The Kinetic-Turbidimetric Assay can be performed in either endotoxin-free glass test tubes or plastic microtiter plates.


An incubating kinetic tube reader or microplate reader and associated software must be used when performing the Kinetic-Turbidimetric Assay (KTA).

### Calculation of endotoxin concentration

During a KTA reaction the turbidity of the test solutions is continually monitored by the Toxinometer®. The time is takes for a sample to reach the predetermined absorbance (Threshold) is known as gelation time. The gelation time is referred to by the Toxinometer® software as (Tg). The software produces a log(x axis)/log(y axis) correlation of the Tg of each standard with its corresponding endotoxin concentration. Below is an example of a representative standard curve.

### **Representative Standard Curve**

Standard	EU/mL	Tg (minutes)
1	0.001	68.8
2	0.01	40.0
3	0.1	23.6
4	1.0	13.2
5	10	8.0
slope	-0.235	
Y intercept	1.13	
correlation	-0.999	



## PYROSTAR<sup>TM</sup> ES-F Series Full Film



PYROSTAR™ ES-F reagents are available with the KTA quantitative ranges noted below. The KTA quantitative range is related to the gel-clot sensitivity.

Gel-Clot Sensitivity (EU/mL)	KTA Quantitative Range (EU/mL)
0.03 - 0.25	0.01 to 10
*0.015	0.001 to 10

<sup>\*</sup>will only test down to 0.001EU/ml when tube reader is used

### **Product Features**

- Endotoxin-specific lysate, avoids false positive results from glucans
- Available in multi-tests vials or single-test vials
- Can be used as either a gel-clot or Kinetic-Turbidimetric Assay (KTA) reagent
- Gel-Clot lysate sensitivities range from 0.015 to 0.25 EU/mL
- KTA assays can be performed in tube reader or microplate reader
- PYROSTAR™ ES-F reagents, when used with a microplate reader, have a KTA quantitative range of 0.01 EU/mL to 10 EU/mL.
- Available with matched control standard endotoxin (CSE)

### PYROSTAR™ ES-F USING MICROPLATE READER

### **Product Features**

- Endotoxin specific testing
- Quantitative detection by Kinetic-Turbidimetric method
- KTA Quantitative Range (EU/mL) 0.01 to 10
- Formulated for no interference from beta-Glucans
- No buffer required to block beta-Glucan interference
- Requires only 50uL of sample compared to 100uL on tube reader

## FUJIFILM PYROSTAR<sup>TM</sup> ES-F Series

### Multi-Test Lysate with matched CSE

PYROSTAR™ ES-F 2.0 mL Test Kit • 4 multi-test vials (2.0 mL) + 1 vial CSE (500 ng/vial)		
Catalog number	Gel-clot Sensitivity EU/mL	Kinetic-Turbidimetric Assay
		Quantitative Range (EU/mL)
WPEK4-20015	0.015	0.001 to 10
WPEK4-20003	0.03	
WPEK4-20006	0.06	0.01 to 10
WPEK4-20125	0.125	
WPEK4-20025	0.25	

PYROSTAR™ ES-F 5.2 mL Test Kit • 4 multi-test vials (5.2 mL) + 1 vial CSE (500 ng/vial)		
Catalog number	Gel-clot Sensitivity EU/mL	Kinetic-Turbidimetric Assay
		Quantitative Range (EU/mL)
WPEK4-50015	0.015	0.001 to 10
WPEK4-50003	0.03	
WPEK4-50006	0.06	0.01 to 10
WPEK4-50125	0.125	
WPEK4-50025	0.25	

### Single- Test Lysate with matched CSE

PYROSTAR™ ES-F Single Test Kit • 25 single-test vials + 1 vial CSE (500 ng/vial)		
Catalog number	Gel-clot Sensitivity EU/mL	Kinetic-Turbidimetric Assay
		Quantitative Range (EU/mL)
WPESK-0015	0.015	0.001 to 10

### **Bulk Lysate**

PYROSTAR™ ES-F 2.0 ml Bulk Kit • 100 multi-test vials (2.0 mL)		
Catalog number	Gel-clot Sensitivity EU/mL	Kinetic-Turbidimetric Assay
		Quantitative Range (EU/mL)
WPEM-20015	0.015	0.001 to 10
WPEM-20003	0.03	
WPEM-20006	0.06	0.01 to 10
WPEM-20125	0.125	
WPEM-20025	0.25	

### FUJIFILM Wako Chemicals U.S.A. Corp. FUJIFILM



### **Bulk Lysate**

PYROSTAR™ ES-F 5.2 ml Bulk Kit • 100 multi-test vials (5.2 mL)		
Catalog number	Gel-clot Sensitivity EU/mL	Kinetic-Turbidimetric Assay
		Quantitative Range (EU/mL)
WPEM-50015	0.015	0.001 to 10
WPEM-50003	0.03	
WPEM-50006	0.06	0.01 to 10
WPEM-50125	0.125	
WPEM-50025	0.25	

### **Our Promise**

As an FDA licensed facility, FUJIFILM Wako Chemicals U.S.A. Corp. - LAL Division is committed to ensuring that our production site and LAL reagents comply with all the rules, regulations, and quality standards set forth by FDA for current Good Manufacturing Practices (cGMP's).

### **Horseshoe Crab Conservation**



FUJIFILM Wako Chemicals U.S.A. Corp. is very much concerned with maintaining the viability of the horseshoe crab population. We are dedicated to following practices that ensure the careful handling and good quality of crabs used for LAL manufacture that both minimize injury and protect this invaluable species. After bleeding, the crabs are returned the next day by our fishermen to the same waters where they were collected. To assist in the collection of data for crab conservation studies, Fujifilm Wako Chemicals U.S.A. Corporation participates in a horseshoe crab tagging and monitoring program coordinated by the U.S. Fish and Wildlife Service.



FUJIFILM Wako Chemicals U.S.A. Corp. - Cape Charles Bleeding Facility: Cape Charles, Virginia 23310

### FUJIFILM Wako Chemicals U.S.A. Corp.

#### **United States of America**

FUJIFILM Wako Chemicals U.S.A. Corp.

1600 Bellwood Road

North Chesterfield, VA 23237 Toll Free: 800-992-9256 Tel: 804-714-1919

Tel: 804-672-4655 Fax: 804-271-7791

Email: wkuspyrostarinfo@fujifilm.com

www.wakopyrostar.com

### Japan

FUJIFILM Wako Pure Chemical Corp. 1-2 Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605 Japan

Telephone: + 81-6-6203-3741 Facsimile: + 81-6-6222-1203 www.wako-chem.co.jp

### Germany

**FUJIFILM Wako Chemicals** Europe GmbH Fuggerstraße 12 41468 Neuss Germany Telephone: + 49-2131-311-0

Facsimile: + 49-2131-311-100 www.wako-chemicals.de

