

The background of the entire page is a scanning electron micrograph (SEM) showing several large, complex, and porous structures. These structures are colored in shades of blue, purple, and pink, set against a dark, almost black background. The structures appear to be highly textured and interconnected, resembling biological cells or synthetic scaffolds. The lighting highlights the intricate details of their surfaces, showing ridges, valleys, and fine filaments.

**Lonza**

Pharma & Biotech

# **PyroCell™ MAT System quick start guide**

# PyroCell™ MAT System

## quick start guide


### Routine use of the PyroCell™ MAT System for *in vitro* detection of pyrogenic contaminants in test samples

The PyroCell™ MAT System, IL-6 comprises cryopreserved PBMC pooled from 4 healthy human donors (pMAT Cells), an optimized MAT Culture Medium Supplement, and a human IL-6 ELISA assay system (00254714). Monocytes, the key cells of innate immunity, respond to pyrogens in a product sample by producing pro-inflammatory cytokines such as IL-6 during a stimulation period. The cell culture supernatant containing the released IL-6 cytokines is then analyzed with the PeliKine compact human IL-6 ELISA assay kit.

Regulatory requirements as outlined for example, by the European Pharmacopeia (Chapter 2.6.30) describe the following options for routine use of the MAT assay: the quantitative test (Method A), the semi-quantitative test (Method B), and the reference lot comparison test (Method C). This quick start guide is describing an example employing Method A or B. The “historical LOD” for the PyroCell™ MAT assay system is 0.02 EU/mL. For more details of the test procedure please refer to the PyroCell™ MAT System User Guide for Routine Use.


#### Day 1: Stimulation of pMAT Cells

- Prepare endotoxin dilutions.
  - Prepare test samples on the cell culture plate.
  - Thaw pMAT Cells, transfer to the plate and incubate for 18-24h.
- 
- Overnight coating of the ELISA plate.




#### Day 2: IL-6 ELISA assay

- Harvest the culture supernatants.
- Detect IL-6 cytokines with an IL-6 ELISA assay.
- Read results in a Microplate absorbance reader.




Calculate the pyrogenicity of the test sample.



# PyroCell™ MAT System

## quick start guide

### Day 1: Stimulation of pMAT Cells



Column 1-3: 5-point standard curve (0.16 – 0.08 – 0.04 – 0.02 – 0.01 EU/mL), and a blank sample.

Column 4-6 & 7-9: test samples (3 dilutions each) with and without control endotoxin spike (0.04 EU/mL).

Column 10-12: test sample or optional NEP control (spiked into complete medium and sample S1-1 and S2-1).

#### Solutions needed

- 1 mL reference standard endotoxin (RSE) dilution, 0.16 EU/mL (by serial dilution).
- 2 mL of endotoxin spike dilution, 0.08 EU/mL (3 samples).
- 2.5 mL of test sample dilutions (up to 3 samples).
- 1 mL non-endotoxin pyrogen (NEP) dilution (optional control, replace the 3rd sample).

#### Points to consider


- **Do not vortex dilutions containing complete medium or pMAT Cells.**
- Prepare all dilutions for cell culture in an aseptic environment (laminar airflow cabinet) at the day of experiment.
- Use pyrogen-free accessories and tubes only. Equilibrate reagents to room temperature before use.
- Avoid bubbles and foam formation.
- The highest dilution on the plate shall not exceed the maximum valid dilution (MVD).
- Final dilution in the plate: 2-fold for all samples.

### Endotoxin dilutions

- Prepare complete medium by adding MAT Culture Medium Supplement to 33 mL IMDM. Mix by inverting the tube 10x.

#### Step 1: Serial dilution

#### Step 2: Spike solution



RSE stock solution in LAL water

50 µL RSE dilution, 950 µL IMDM

20 µL RSE dilution, 980 µL complete medium


80 µL RSE dilution, 920 µL complete medium

80 µL RSE solution 1920 µL complete medium

# PyroCell™ MAT System


## quick start guide

### Prepare the cell culture plate




#### Step 3: Transfer complete medium

- Add 100 µL of complete medium into every yellow well.



#### Step 4: Prepare the endotoxin standard curve

- Add 200 µL of the 0.16 EU/mL RSE dilution into wells A1-D1.
- Transfer 100 µL from A1-D1 to A2-D2. Mix by pipetting up and down 10x.
- Repeat transfer from A2-D2 to A3-D3, then to E1-H1 and E2-H2 as illustrated.
- Discard 100 µL from E2-H2.




#### Step 5: Prepare the sample dilution

- Add 200 µL of the sample S1 into column 4.
- Transfer 100 µL from column 4 to 5, and column 5 to 6. Mix each time by pipetting up and down 10x.
- Discard 100 µL from column 6.
- Repeat 2-fold dilution for other test samples, e.g. sample S2 (columns 7-9).
- [Optional for NEP] Add 100 µL of sample S1 into E10-H10 and sample S2 into E11-H11.

# PyroCell™ MAT System

## quick start guide



### Step 6: Add endotoxin spike and optional, NEP control

- Add 50  $\mu$ L of complete medium into yellow wells.
- Add 50  $\mu$ L of endotoxin spike solution (0.08 EU/mL) to the grey wells.
- Add 50  $\mu$ L of NEP to column 10 and 11 (striped wells).

## Thaw pMAT Cells, transfer to plate and incubate




### Step 1

Thaw a vial of pMAT Cells at 37°C until observing a small remaining clump of ice.




### Step 2

Immediately transfer thawed cells into a 50 mL tube and slowly add 5 mL of complete medium while swirling gently. **Do not vortex!**



### Step 3

Fill cell suspension into a reservoir and mix by pipetting up and down.



### Step 4

Transfer 50  $\mu$ L pMAT Cells to each well of columns 1-11; incubate the plate at 37°C for 18-24h.

# PyroCell™ MAT System

## quick start guide


### Day 2: IL-6 ELISA assay

#### Points to consider

- During wash steps make sure to empty all wells completely after the last wash step.
- Equilibrate all buffers and solutions to room temperature prior to use (exception: anti-IL 6 antibody and streptavidin-HRP reagent).

#### Step1: Overnight coating of the ELISA plate (day 1)

- Dissolve content of one coating buffer capsule in 100 mL distilled water, incubate 5 minutes.
- Mix 12 mL of coating buffer with 120 µL coating antibody.
- Add 100 µL coating buffer to each wells of the microtiter plate and incubate over night at room temperature.




Step 2



#### Harvesting the cell culture supernatants

- Transfer supernatants from the cell culture plate (day1) into a fresh plate (plate 1).
- Dissolve 15 mL 5-fold concentrated HPE buffer in 60 mL distilled water.
- Prepare a dilution of each supernatant (e.g. a 1:5 dilution) by adding 120 µL of 5-fold concentrated HPE buffer into each well of a fresh plate (plate 2). Transfer 30 µL supernatant from plate 1 to the corresponding well of plate 2. Mix by pipetting up and down.






Step 3



|            |   |  |
|------------|---|--|
| ELISA step | Wash  | Blocking of unspecific binding sites                   |
| Reagent    | PBS Buffer (1 PBS tablet in 200 mL distilled water) | Blocking Buffer (500 µL Blocking Reagent in 25 mL PBS) |
| Use        | 300 µL/ well<br>4 times                             | 200 µL/ well<br>1h incubation at 18-25°C               |

# PyroCell™ MAT System

## quick start guide

|        |   |  |  |   |  |
|--------|---|--|--|---|--|
| Step 4 |    | <p><b>ELISA step</b></p> <p><b>Reagent</b></p> <p><b>Process</b></p>   | <p>Wash</p> <p>Wash Buffer (50 mL concentrate in 950 mL distilled water)</p> <p>300 µL/ well<br/>5 times</p> | → | <p>Binding of Cytokines from supernatant (Optional, add IL-6 Standard Dilutions)</p> <p>Supernatant dilutions (plate 2)</p> <p>100 µL/ well<br/>1h incubation at 18-25°C</p> |
| Step 5 |    | <p><b>ELISA step</b></p> <p><b>Reagent</b></p> <p><b>Process</b></p>   | <p>Wash</p> <p>Wash Buffer (from Step 3)</p> <p>300 µL/ well<br/>5 times</p>                                 | → | <p>Binding of anti-cytokine antibody</p> <p>Biotinylated anti-IL-6 antibody (120µL antibody + 12 mL HPE)</p> <p>100 µL/ well<br/>1h incubation at 18-25°C</p>                |
| Step 6 |  | <p><b>ELISA step</b></p> <p><b>Reagent</b></p> <p><b>Process</b></p>   | <p>Wash</p> <p>Wash Buffer (from Step 3)</p> <p>300 µL/ well<br/>5 times</p>                                 | → | <p>Binding of HRP conjugate</p> <p>Streptavidin-HRP onjugated antibody (1.2 µL antibody + 12 mL HPE)</p> <p>100 µL/ well<br/>30 min incubation at 18-25°C</p>                |
| Step 7 |  | <p><b>ELISA step</b></p> <p><b>Reagent</b></p> <p><b>Process</b></p>   | <p>Wash</p> <p>Wash Buffer (from step 3)</p> <p>300 µL/ well<br/>5 times</p>                                 | → | <p>Color reaction</p> <p>TMB substrate, 100 µL</p> <p>Up to 10 min incubation (dark) at 18-25°C<br/>Add Stop Solution, 100 µL</p>  |
| Step 8 |  | <p><b>Read results</b></p> <ul style="list-style-type: none"> <li>Determine the absorptions (OD) at 450 nm and 540-590 nm in a plate reader (e.g. ELx808 Reader) within 30 minutes.</li> <li>Analyze results.</li> </ul> |  |   |  |

## MAT reagents and pyrogen-free materials

| Product name                                      | Cat. No. |
|---|----------|
| PyroCell™ MAT System, IL-6                        | 00254714 |
| PyroCell™ MAT System, IL-1 beta                   | 00254715 |
| Iscove's Modified Dulbecco's Medium (IMDM)        | 12-722F  |
| Reference Standard Endotoxin (RSE – USP)          | E700     |
| LAL Reagent Water (e.g. 100 mL)                   | W50-100  |
| Sterile, pyrogen-free glass tubes                 | N207     |
| Reagent Reservoir                                 | 00190039 |
| LAL Reagent Grade™ Multi-well Plates              | 25-340   |
| Pyrogen-free Eppendorf® Biopur Tips, 2-200 µL     | 25-415   |
| Pyrogen-free Eppendorf® Biopur Tips, 50-1000 µL   | 25-417   |
| ELx808™ Reader with filters (450 nm & 540-590 nm) | 25-315S  |
| (optional) Non-endotoxin pyrogens                 |          |

### What you also need

- 50 mL tubes (cell-culture grade, endotoxin-free)
- Bottles and test tubes (buffer preparation)
- Volumetric & serological pipette (as needed)
- Adjustable pipettes, multichannel pipette,
- pipette tips (standard)
- 96-well round bottom plates (standard)
- Laminar airflow cabinet (aseptic)
- CO<sub>2</sub> cell incubator (37°C, 5% CO<sub>2</sub>)
- Water bath (37°C)
- Automated microplate washer and microplate shaker

### Supporting Documents

- Certificate of Analysis (CoA), [www.lonza.com/coa](http://www.lonza.com/coa)
- PyroCell™ MAT System User Guide for Routine Use
- PyroCell™ MAT Analytics Guide

### Contact information

#### North America

Customer Service: +1 800 638 8174 (toll free)  
Fax: +1 301 845 8338  
[order.us@lonza.com](mailto:order.us@lonza.com)  
Scientific Support: +1 800 521 0390 (toll free)  
[scientific.support@lonza.com](mailto:scientific.support@lonza.com)

#### Europe

Customer Service: +32 87 321 611  
[order.europe@lonza.com](mailto:order.europe@lonza.com)  
Scientific Support: +49 221 99199 400  
[scientific.support.eu@lonza.com](mailto:scientific.support.eu@lonza.com)

#### International

Contact your local Lonza distributor  
Customer Service: +1 301 898 7025  
[scientific.support@lonza.com](mailto:scientific.support@lonza.com)

Lonza Walkersville, Inc. – Walkersville, MD 21793

For research use only. Not for use in diagnostic procedures.

All trademarks belong to Lonza or its affiliates or to their respective third party owners. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights. The user bears the sole responsibility for determining the existence of any such third party rights, as well as obtaining any necessary licenses.

@2020 Lonza. All rights reserved.  
RT-MN026 05/20

[www.lonza.com/mat](http://www.lonza.com/mat)  
[www.bioscience.lonza.com](http://www.bioscience.lonza.com)