

# Endotoxin detection – which method is best for me?

# Endotoxin detection methods

## LAL Assays

- Gel clot
- Semi-quantitative gel clot
- *Endpoint chromogenic*
- *Endpoint turbidimetric*
- Kinetic chromogenic
- Kinetic turbidimetric

## Recombinant Factor C Assays

- Endpoint fluorescent



# Which test is best?



Each test has its own place



The choice of test for a laboratory depends on:

- ➔ Result required
  - Quantitative/Semi-quantitative
  - Qualitative

---

- ➔ Available equipment

---

- ➔ Sample to be tested



# Lonza PYROGENT® Gel Clot LAL Assay


- PYROGENT® Kits
  - lysate only (matched, lyophilized control standard endotoxin (CSE) can be purchased separately)
- PYROGENT® Plus Kits
  - lysate and matched, lyophilized CSE
- Multiple sensitivities and vial sizes available

Vial size	Available sensitivities (EU/mL)
50 tests/vial	0.03, 0.06, 0.125, 0.25
16 tests/vial	0.06, 0.125


- Multiple kit configurations – standard sizes and bulk



# Gel clot assay



# Performing the gel clot assay



# Reading the gel clot assay

- Invert tube 180° in a smooth motion
- A positive reaction is characterized by the formation of a firm gel that remains intact
- A negative reaction is characterized by the absence of a solid gel after inversion



# PYROGENT<sup>®</sup> Gel Clot LAL Assay

**LONZA**

Pharma & Biotech



## Advantages

- A simple test to perform
- Low start-up costs
- No readers/software necessary



## Disadvantages

- Results are subjective
- Gel formation is prone to interference
- Labor intensive
- Qualitative  
Semi-quantitative with many dilutions





# Gel clot assays



The gel clot assay can be performed in two ways:

Limit test

Semi-quantitative

# Gel clot assays

## Limit Test

- Yes/No answer to a specific kit sensitivity (more than/less than labeled lysate sensitivity)

---

- Positive control, negative control, sample and sample positive product control (PPC)

## Semi-quantitative

- Standard series of endotoxin

---

- Dilution series of sample

---

- PPC most concentrated sample

---

- Negative control

---

- Sample endotoxin load range = Dilution factor of last positive (clotting) sample x kit sensitivity



# Kinetic LAL assays

01

The most recently developed LAL assays

02

Most objective of all the LAL assays

03

Two forms of kinetic LAL assays are available:

- Kinetic turbidimetric
- Kinetic chromogenic




# Lonza kinetic LAL assays


- Kinetic turbidimetric – Lonza PYROGENT®-5000 Assay
- Kinetic chromogenic – Lonza Kinetic-QCL® Assay
- Both kinetic LAL assay types are quantitative and can be conducted using an eight-channel incubating plate reader (ELx808™ Reader) with supporting software (Lonza WinKQCL® Software), provided the correct filters are used:
  - PYROGENT®-5000 Assay – 340 nm filter
  - Kinetic-QCL® Assay – 405 nm filter




# Kinetic assays



# Performing the kinetic assay



# How kinetic assays work



- The optical density in each well is read at the start of the assay

---

- All results from the first reading are set to a baseline level to account for color, opacity and more

---

- The reader measures the optical density at fixed time intervals looking for a delta mOD of 30 for the kinetic turbidimetric assay and a delta mOD of 200 for the kinetic chromogenic assay


---

- The 'end of assay' mark is reached when the lowest standard has reached the required mOD change

---

- The time taken for each well to reach the required mOD change is known as the reaction time

# Kinetic LAL reaction





# The kinetic assay



The mean reaction time for each of the endotoxin standards is plotted against endotoxin content using linear regression or PowerCurve™ to create a standard curve




The non-linear nature of the data means that at a minimum, the data must be subject to a log-log transformation to plot the standard curve



Mean reaction times for the sample and sample PPCs are used to calculate interpolated endotoxin content from this standard curve and reported in the relevant units – EU/ml or EU/mg – according to the endotoxin limit entered into the product details



# Reaction time versus endotoxin concentration



A standard curve for the kinetic turbidimetric assay.

Note the log-log ratio of the standards and also the curve slopes downward – as more endotoxin equals a faster assay reaction time.

# PYROGENT<sup>®</sup>-5000 Kinetic Turbidimetric LAL Assay



## Advantages

- Quantitative over a large range – 0.01 to 100 EU/ml
- Good for testing water and simple products
- Lower reagent cost than kinetic chromogenic assay
- Can cope with yellow-colored products which absorb in the 405 nm range (i.e. unsuitable for chromogenic testing)



## Disadvantages

- Less sensitive than the Kinetic-QCL<sup>®</sup> Assay
- Similar interference problems as in the gel clot method and not suitable for turbid or viscous products
- Prone to problems from bubbles
- Requires a reader, computer and software



# Kinetic-QCL<sup>®</sup> Kinetic Chromogenic LAL Assay



## Advantages

- Less subject to interference and copes with turbid or more viscous products

---

- Lonza's most sensitive LAL assay

---

- Quantitative over a large range – 0.005 to 50 EU/ml



## Disadvantages

- Requires a reader, computer and software

---

- Not suitable for colored products which absorb in the 405 nm range




# A sustainable alternative to LAL – PyroGene® Assay

- PyroGene® Recombinant Factor C Assay is an animal-free, sustainable alternative to LAL-based assays
- This assay has several major advantages over conventional LAL assays
- The future of endotoxin testing



# Recombinant factor C endotoxin detection assay

## *Recombinant Factor C*



# PyroGene<sup>®</sup> Assay overview

- Single step quantitative endpoint fluorescent assay
- Assay takes one hour at 37°C
- Uses liquid reagents
- Performed in standard 96-well microplates
- The PyroWave<sup>®</sup> XM Reader is controlled by WinKQCL<sup>®</sup> Software




# PyroGene® Assay features

- Consistent assay performance
- Security of supply
  - No animal utilization
- Specificity for endotoxin
  - Elimination of glucan reaction pathway removing false positive issue seen with LAL
- Comparability to LAL-based methods in terms of performance and pricing






# Performing the PyroGene<sup>®</sup> Assay




Reconstitute CSE as per package insert and prepare a 4-point standard curve 0.005 to 5 EU/ml




Prepare standards, samples and controls and pipette into 96-well plate



Pre-incubate the plate for 10 minutes at 37°C



Prepare working reagent according to package insert (ratio 5:4:1, substrate:buffer:enzyme) and add to each well



Start assay

# Reading the PyroGene<sup>®</sup> Assay

01

WinKQCL<sup>®</sup> Software will take the first read for all wells and store them

02

The plate is then incubated for 60 minutes at 37°C

03

At the end of the assay the plate is read again, and the readings stored



# PyroGene<sup>®</sup> Assay – calculations

WinKQCL<sup>®</sup> Software  
then calculates:


- ➔  $\Delta\text{RFU} = \text{RFU}_{60 \text{ min}} - \text{RFU}_{\text{time 0}}$
- ➔ Averages  $\Delta\text{RFU}$  for BLANKs
- ➔ Subtracts mean  $\Delta\text{RFU}_{\text{blank}}$  from all  $\Delta\text{RFUs}$

WinKQCL<sup>®</sup> Software  
then constructs a  
log-log plot of  $\Delta\text{RFU}$   
versus endotoxin  
concentration

Unknowns (samples  
and sample PPCs) are  
determined using this  
standard curve






# PyroGene<sup>®</sup> rFC Assay standard curve



A standard curve for the PyroGene<sup>®</sup> rFC Assay.

The standard curve has a low RFU for the 0.005 EU/ml standard and a high RFU for the 5 EU/ml standard.

# Endotoxin assay overview



- PYROGENT®-5000 Kinetic Turbidimetric LAL Assay
- Kinetic-QCL® Kinetic Chromogenic LAL Assay
- PyroGene® Recombinant Factor C Assay

- 96-well plates
- Pipette tips
- Reagent reservoirs
- Pyrogen-free dilution tubes

- Fully integrated solution for endotoxin detection and analysis
- Built with components to meet the latest data integrity requirements

# Summary of endotoxin detection methods

- LAL Assays
  - Gel clot
  - Semi-quantitative gel clot
  - *Endpoint chromogenic*
  - *Endpoint turbidimetric*
  - Kinetic chromogenic
  - Kinetic turbidimetric
- Recombinant Factor C Assays
  - Endpoint fluorescent



# Choice of assay is dependent on...

	Gel clot	Kinetic turbidimetric	Kinetic chromogenic	rFC assay
Workload	● ● ●	●	●	●
Nature of the samples to be tested	Water, in-process and final release testing. Also suitable for plant-based material and acidic/basic samples	Water samples, simple products and yellow-colored products which absorb in the 405 nm range	Biological products (i.e. vaccines and antibiotics) and turbid or viscous samples	Water, in-process, final release testing. Also suitable for plant-based material
Available equipment or equipment budget	Dry heat block or water bath	Incubating absorbance reader	Incubating absorbance reader	Incubating fluorescence reader
Need for quantitation	Qualitative (Yes/No answer)	Quantitative (Results calculated from a standard curve)	Quantitative (Results calculated from a standard curve)	Quantitative (Results calculated from a standard curve)

# Thank You

ELx808 is a trademark of BioTek Instruments, Inc.

All trademarks belong to Lonza, registered in the USA, EU or CH or to third party owners and used only for informational purposes. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. For more details: [www.lonza.com/legal](http://www.lonza.com/legal). © 2021 Lonza. RT-PP025 06/2021