

Transfection of Human Induced Pluripotent Stem Cells using Nucleofector™ Technology

Technical Reference Guide

Introduction

For studying differentiation pathways in human induced pluripotent stem cells (iPSCs) or modifying such pathways, it is important to find a transfection method which is efficient while maintaining pluripotency. The non-viral Nucleofector™ Technology has proven to be well suited for both iPSC generation¹⁻⁴ as well as iPSC transfection, e.g. for genome editing applications⁵⁻¹¹. Here we present general guidelines for establishing Nucleofection™ Conditions for transfecting iPSCs. For more specific recommendations about genome editing please visit www.lonza.com/genome-editing.

1. Basic Protocol Concept

Similar to human embryonic stem cell (ESC) clones, generated iPSC clones may differ clone to clone regarding optimal Nucleofection™ Conditions. Thus, providing a ready-to-use Optimized Protocol with pre-defined conditions (i.e. one solution + one program) which are guaranteed to work with any clone is not possible. For each iPSC clone, it is recommended to first determine the optimal conditions using our pmaxGFP™ Control Vector following the „Basic Human Stem Cell Protocol“ for the respective Nucleofection™ Platform*. A basic protocol describes how to easily define optimal Nucleofection™ Conditions by testing a pre-defined matrix of Nucleofector™ Programs and Solutions.

*For downloading the platform specific detailed protocol, please refer to www.lonza.com/protocols and filter for “stem cells”.

2. Optimization Process

2.1 First Optimization Round

For 4D-Nucleofector™ X Unit:** When using the 4D-Nucleofector™ System the initial optimization experiment is comprised of 32 reactions. Seven different Nucleofector™ Programs plus a no program control are tested in duplicate with two 4D-Nucleofector™ Solutions (Table 1).

** The same applies for the 96-well Shuttle™ Add-on

Table 1. Optimization scheme for 4D-Nucleofector™ System

P3 Solution (Strip 1, n=2)		P4 Solution (Strip 2, n=2)	
CA-137	CA-137	CA-137	CA-137
CB-150	CB-150	CB-150	CB-150
CD-118	CD-118	CD-118	CD-118
CE-118	CE-118	CE-118	CE-118
CM-113	CM-113	CM-113	CM-113
DC-100	DC-100	DC-100	DC-100
DN-100	DN-100	DN-100	DN-100
No program	No program	No program	No program

For Nucleofector™ II/2b Device: When using the Nucleofector™ II/2b Device, the initial optimization experiment is comprised of 12 reactions. Five different Nucleofector™ Programs plus a no program control are tested in combination with two Nucleofector™ Solutions (Table 2).

Table 2. Optimization scheme for Nucleofector™ II/2b Device

HSC ₁ Solution	HSC ₂ Solution
No program	No program
A-012	A-012
A-013	A-013
A-023	A-023
A-027	A-027
B-016	B-016

The Nucleofector™ Program and Solution which turned out to be the most appropriate Nucleofection™ Condition should be used for all subsequent transfections.

2.2 Fine Tuning of Conditions

If results of the first optimization round are not fully satisfactory, our Scientific Support Team can provide fine-tuning suggestions.

Depending on the clone, transfection efficiencies up to 60% and viabilities between 50-75% can be achieved. For more details on results, you may check our cell transfection database (www.lonza.com/celldata-base) or the various publications⁵⁻¹¹.

3. Factors Influencing Nucleofection™ Process Results

3.1 Cell Handling Prior to Transfection

Single cell suspension: We recommend transfecting the cells in single cell suspension. Transfection of clumps leads to lower transfection efficiency and less reproducibility¹². If single cell suspension passaging is not established, please perform some pre-experiments by testing Accutase^{13, 14} and Trypsin¹⁵ for detachment. Cultivate the cells afterwards and analyze which method led to highest viability and lowest differentiation.

Detachment of stem cells with feeder cell culture: There are three possibilities to remove feeder cells from your stem cell culture prior to Nucleofection™ Process:

- If stem cells are usually cultured on feeder cells, passage them once in matrix coated plates to remove the feeder cells (for details see reference 2). Then proceed as with feeder-free cultured cells.
- Cultivate the cells on feeder cells until the day of the experiment. Detach the stem cells with Collagenase. Then dissociate the clumps with Accutase into a single cell suspension.
- Cultivate the cells on feeder cells until the day of the experiment. Detach all cells with Accutase. Incubate the cells on an uncoated cell culture flask for 1 hour in a humidified 37°C/5% CO₂ incubator. The feeder cells will attach and the stem cells will stay in suspension. Harvest the cells in suspension.

Detachment of stem cells with feeder-free cultures: Prior to Nucleofection™ Process detach the cells from the matrix-coated plates by incubation with Accutase for 5 minutes at 37°C. Dissociate the cells into a single cell suspension by pipetting the suspension carefully up and down 4 – 6 times. Add medium to stop the Accutase reaction.

Apoptosis inhibitors: The use of apoptosis inhibitors like ROCK inhibitor¹⁶ and neurotrophins¹⁷ prior to and post-transfection have been reported to increase viability of human ESCs. For iPSCs, using apoptosis inhibitors seems to be advantageous as well^{5,11}.

3.2 Cell Handling Post-Transfection

Plating Density: The plating density post-transfection is a critical aspect for the viability of human pluripotent stem cells. Our experience shows that higher densities lead to better viability of the cells. Therefore, we recommend plating human pluripotent stem cells at densities from 4 - 6.5 x 10⁵ cells per cm². These plating cell numbers and volumes produce optimal Nucleofection™ Results in most cases. However, an extended range of cell numbers might be tested depending on your specific needs. Because of the high plating density, a passaging step 48 hours post transfection using collagenase or dispase might be required.

Culture on matrix: When cells are cultured on matrix-coated plates the coating should be fresh. Storage of coated plates for more than 7 days leads to reduced attachment of the cells post-transfection. After cells have been transferred from the Nucleofection™ Vessel to the culture plates, the plates should be centrifuged for 3 minutes at 70xg (at room temperature).

References

- 1) Chen *et al* (2013) PLoS ONE 8(10): e75682
- 2) Goyal *et al* (2013) PLoS ONE 8(12): e82838
- 3) Chou BK *et al* (2011) Cell Research 21:518-529
- 4) Hu K *et al* (2011) Blood 117(14): e109-e119
- 5) Chang C.J *et al* (2012) Blood 120:3906 - 3911
- 6) Ou W *et al* (2013) PLoS ONE 8(11): e81131
- 7) Yang L *et al* (2013) Nucleic Acids Res 41:9049 – 9061
- 8) Zhu F *et al* (2013) Nucleic Acids Res 10: 1093
- 9) Zou J *et al* (2009) Stem Cells 5(1):97-110
- 10) Zou J *et al* (2011) Blood 117(21):5561-5572
- 11) Zou J *et al* (2011) Blood 118(17):4599-4608
- 12) Hohenstein KA *et al* (2008) Stem Cells 26(6):1436-43
- 13) Siemen H *et al* (2005) Stem Cells Dev 14:378-383
- 14) Bajpai R *et al* (2008) Mol Reprod Dev 75(5):818-27
- 15) Ellerström C *et al* (2007) Stem Cells 25:1690-1696
- 16) Watanabe K *et al* (2007) Nat Biotech 25 (6): 681-686
- 17) Pyle AD *et al* (2006) Nat Biotech 24(3): 344-350

Contact Information

North America

Customer Service: +1 800 638 8174 (toll free)
order.us@lonza.com
Scientific Support: +1 800 521 0390 (toll free)
scientific.support@lonza.com

Europe

Customer Service: +32 87 321 611
order.europe@lonza.com
Scientific Support: +32 87 321 611
scientific.support.eu@lonza.com

International

Contact your local Lonza distributor
Customer Service: +1 301 898 7025
Fax: +1 301 845 8291
scientific.support@lonza.com

International Offices

Australia	+ 61 3 9550 0883
Belgium	+ 32 87 321 611
Brazil	+ 55 11 2069 8800
France	0800 91 19 81 (toll free)
Germany	0800 182 52 87 (toll free)
India	+91 40 4123 4000
Japan	+ 81 3 6264 0660
Luxemburg	+32 87 321 611
Singapore	+ 65 6521 4379
The Netherlands	0800 022 4525 (toll free)
United Kingdom	0808 234 97 88 (toll free)

Lonza Cologne GmbH – 50829 Cologne Germany

For research use only. Not for use in diagnostic procedures.
The Nucleofector™ Technology is covered by patent and/or patent pending rights owned by the Lonza Group Ltd or its affiliates.
Unless otherwise noted, all trademarks herein are marks of the Lonza Group Ltd or its affiliates. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights.

© Copyright 2014, Lonza Cologne GmbH. All rights reserved.
CD-DS020 12/14
