


Amaxa® 4D-Nucleofector® Protocol for Immature Mouse Dendritic Cells For 4D-Nucleofector® X Unit — Transfection in suspension

Progenitor cells derived from bone marrow; non-adherent or loosely adherent cells of irregular shape with typical protrusions ("dendrites") of variable shape and length

Example for Nucleofection® of immature mouse dendritic cells

Average transfection efficiency of immature mouse dendritic cells 24 hours post Nucleofection®. 5×10^4 immature mouse dendritic cells were transfected with program DK-100 (high efficiency) or DJ-100 (high viability) and $0.4 \mu g$ of pmaxGFP® Vector in $20 \mu l$ Nucleocuvette® Strips. 24 hours post Nucleofection® cells were analyzed on a FACSCalibur™ [Becton Dickinson]. Cell viability was approximately 74% for DK-100 and 98% for DJ-100 after 24 hours [% PI-negative compared to untreated].



Product Description

Recommended Kit(s) – P4 Primary Cell 4D-Nucleofector® X Kit

Cat. No.	V4XP-4012	V4XP-4024	V4XP-4032
Transfection volume	100 μΙ	100 μΙ	20 μΙ
Size [reaction]	2 x 6	24	2 x 16
Nucleofector® Solution	2 x 0.675 ml (0.492 ml + 27% overfill)	2.25 ml (1.968 ml + 13% overfill)	0.675 ml (0.525 ml + 22% overfill)
Supplement	2 x 0.15 ml (0.108 ml + 27% overfill)	0.5 ml (0.432 ml + 13% overfill)	0.15 ml (0.115 ml + 22% overfill)
pmaxGFP $^{\circ}$ Vector (0.5 μ g/ μ l in 10 mM Tris pH 8.0)	30 µg	30 µg	
pmaxGFP $^{\odot}$ Vector (0.2 μ g/ μ l in 10 mM Tris pH 8.0)			_45 μg
Single Nucleocuvette® (100 μl)	12	24	
16-well Nucleocuvette® Strips (20 µl)	-	-	2

Storage and stability

Store Nucleofector® Solution, Supplement and pmaxGFP® Vector at 4° C. For long-term storage, pmaxGFP® Vector is ideally stored at -20°C. The expiration date is printed on the solution box. Once the Nucleofector® Supplement is added to the Nucleofector® Solution, it is stable for three months at 4° C.

Note

4D-Nucleofector® Solutions can only be used with Nucleovettes® (conductive polymer cuvettes), i.e. in the 4D-Nucleofector® System and the 96-well Shuttle® Device. They are not compatible with the Nucleofector® II/2b Device.

Required Material

Note

Please make sure that the entire supplement is added to the Nucleofector $^{\circ}$ Solution. The ratio of Nucleofector $^{\circ}$ Solution to supplement is 4.5 : 1 (see table 1)

- 4D-Nucleofector® System (4D-Nucleofector® Core Unit and 4D-Nucleofector® X Unit)
- Supplemented 4D-Nucleofector® Solution at room temperature
- Supplied 100 μl single Nucleocuvette® or 20 μl 16-well Nucleocuvette®
 Strips
- Compatible tips for 20 µl Nucleocuvette® Strips: epT.I.P.S. [US/CDN: Eppendorf North America, Cat. No. 2491.431, Rest of World: Eppendorf AG, Cat. No. 0030073.266], Matrix TallTips® [Matrix Technologies Corp., Cat. No. 7281] or LTS Tips [Rainin Instruments, LLC, Cat. No. SR-L10F, SR/SS-L250S, SR/SS-L300S]. Before using other types of pipette tips, please ensure they reach the bottom of the Nucleocuvette® Wells without getting stuck
- Supplied pmaxGFP® Vector
- Substrate of interest, highly purified, preferably by using endotoxin-free kits; A260:A280 ratio should be at least 1.8

Note

As contamination of the DNA with low molecular weight compounds may reduce cell viability, we strongly recommend performing an additional purification step after using a plasmid purification kit. Therefore precipitate the purified DNA twice using 20% PEG/2.5 M NaCl (final concentration).

- Cell culture plates of your choice
- Culture medium I: RPMI 1640 [Lonza; Cat. No. 12-167F] supplemented with 10% calf serum (FCS), 100 μg/ml streptomycin, 100 U/ml Penicilin, 2 mM glutamine and 2000 U/ml GM-CSF (BD Pharmingen; Cat. No.:554586)
- Culture medium II: Culture medium I without GM-CSF
- Prewarm appropriate volume of culture medium I to 37°C (see table 2)
- Appropriate number of cells/sample (see table 2)

1. Pre Nucleofection®

Preparation of cells and cell culture

Preparation of bones

1.1 Carefully remove the femurs and tibia of freshly prepared mouse hind legs using forceps and scissors. Cut off both ends of the bones

Isolation of dendritic cells

- 1.2 Use a 27G needle mounted to a 5 ml syringe to gently flush the bone marrow into a petri dish. Use 2 3 ml culture medium II per bone
- 1.3 Count the viable cells
- 1.4 Spin down cells at 300 xg for 10 minutes at RT and discard the supernatant
- 1.5 Resuspend the cell pellet in culture medium I to reach a cell density of 1 x 10^6 cells/ml
- 1.6 Transfer the cells into 24-well plates (1 ml/well) and incubate them in a 37°C incubator with a 5% CO₃ atmosphere

Note

To yield a high number of functional dendritic cells it is necessary to maintain a sufficient level of GM-CSF. Fresh culture medium I (containing GM-CSF) should be added every second day.

- 1.7 On day 2, carefully remove 700 μ l culture medium I from each well and replace it by fresh culture medium I, to maintain an appropriate GM-CSF concentration
- 1.8 On day 3, remove and discard the culture medium I completely. Wash the cells carefully with 500 µl per well using culture medium II to remove residual non adherent cells
- 1.9 Add 1 ml fresh culture medium I per well
- 1.10 Incubate the cells at 37°C in an incubator with 5% CO₂ atmosphere
- 1.11 On day 6, harvest the dendritic cells by collecting non adherent cells and loosely adherent cells. To release loosely adherent cells, wash off the cells thoroughly by pipetting with culture medium I. Discard the adherent cells
- 1.12 Continue with Nucleofection® as described in chapter 2

2. Nucleofection®

For Nucleofection® Sample contents and recommended Nucleofector® Program, please refer to Table 3.

- 2.1 Please make sure that the entire supplement is added to the Nucleofector® Solution
- 2.2 Start 4D-Nucleofector® System and create or upload experimental parameter file (for details see device manual)
- 2.3 Select/Check for the appropriate Nucleofector® Program (see table 3)
- 2.4 Prepare cell culture plates by filling appropriate number of wells with desired volume of recommended culture media (see table 4) and pre-incubate/equilibrate plates in a humidified 37°C/5% CO₂ incubator
- 2.5 Pre-warm an aliquot of culture medium to 37°C (see table 4)
- 2.6 Prepare plasmid DNA or pmaxGFP® Vector or siRNA (see table 3)
- 2.7 Count an aliquot of the cells and determine cell density

- 2.8 Centrifuge the required number of cells (see table 3) at 300xg for 10 minutes at room temperature. Remove supernatant completely
- 2.9 Resuspend the cell pellet carefully in room temperature 4D-Nucleofector® Solution (see table 3)
- 2.10 Prepare mastermixes by dividing cell suspension according to number of substrates
- 2.11 Add required amount of substrates to each aliquot (max. 10% of final sample volume)
- 2.12 Transfer mastermixes into the Nucleocuvette® Vessels

Note

As leaving cells in Nucleofector® Solution for extended periods of time may lead to reduced transfection efficiency and viability it is important to work as quickly as possible. Avoid air bubbles while pipetting.

- 2.13 Gently tap the Nucleocuvette® Vessels to make sure the sample covers the bottom of the cuvette
- 2.14 Place Nucleocuvette® Vessel with closed lid into the retainer of the 4D-Nucleofector® X Unit. Check for proper orientation of the Nucleocuvette® Vessel
- 2.15 Start Nucleofection® Process by pressing the "Start" on the display of the 4D-Nucleofector® Core Unit (for details, please refer to the device manual)
- 2.16 After run completion, carefully remove the Nucleocuvette® Vessel from the retainer
- 2.17 Resuspend cells with pre-warmed medium (for recommended volumes see table 5). Mix cells by gently pipetting up and down two to three times. When working with the 100 µl Nucleocuvette® use the supplied pipettes and avoid repeated aspiration of the sample
- 2.18 Plate desired amount of cells in culture system of your choice (for recommended volumes see table 5)

3. Post Nucleofection®

3.1 Incubate the cells in humidified $37^{\circ}\text{C}/5\%$ CO₂ incubator until analysis. Gene expression or down regulation, respectively, is often detectable after only 4-8 hours

Additional Information

For an up-to-date list of all Nucleofector® References, please refer to: www.lonza.com/nucleofection-citations

For more technical assistance, contact our Scientific Support Team:

USA /Canada

Phone: 800 521 0390 (toll-free)

Fax: 301 845 8338

E-mail: scientific.support@lonza.com

Europe and Rest of World

Phone: +49 221 99199 400 Fax: +49 221 99199 499

E-mail: scientific.support.eu@lonza.com

References

1. Inaba K. et al, [1992] J. Exp. Med. 176: 1693-1702

Lonza Cologne AG 50829 Cologne, Germany

Please note that the Amaxa® Nucleofector® Technology is not intended to be used for diagnostic purposes or for testing or treatment in humans.

The Nucleofector® Technology, comprising Nucleofection® Process, Nucleofector® Device, Nucleofector® Solutions, Nucleofector® 96-well Shuttle® System and 96-well Nucleocuvette® plates and modules is covered by patent and/or patent-pending rights owned by Lonza Cologne AG.

Amaxa, Nucleofector, Nucleofection, 4D-Nucleofector, Nucleocuvette and maxGFP are registered trademarks of the Lonza Cologne AG in Germany and/or U.S. and/or other countries.

Other product and company names mentioned herein are the trademarks of their respective owners.

This kit contains a proprietary nucleic acid coding for a proprietary copepod fluorescent protein intended to be used as a positive control with this Lonza product only. Any use of the proprietary nucleic acid or protein other than as a positive control with this Lonza product is strictly prohibited. USE IN ANY OTHER APPLICATION REQUIRES A LICENSE FROM EVROGEN. To obtain such a license, please contact Evrogen at license@evrogen.com. The CMV promoter is covered under U.S. Patents 5,168,062 and 5,385,839 and its use is permitted for research purposes only. Any other use of the CMV promoter requires a license from the University of Iowa Research Foundation, 214 Technology Innovation Center, Iowa City, IA 52242.

The use of this product in conjunction with materials or methods of third parties may require a license by a third party. User shall be fully responsible for determining whether and from which third party it requires such license and for the obtainment of such license.

No statement is intended or should be construed as a recommendation to infringe any existing patent.

© Copyright 2010, Lonza Cologne AG. All rights reserved.

Table 1: Volumes required for a single reaction

	100 μl Single Nucleocuvette®	20 μl Nucleocuvette® Strip
Volume of Nucleofector® Solution	82 µl	16.4 µl
Volume of Supplement	18 µl	3.6 µl

Table 2: Required amounts of cells and media for Nucleofection®

	100 µl Single Nucleocuvette®	20 μl Nucleocuvette® Strip
Culture medium per sample post Nucleofection® (for transfer and culture)	0.8 ml	190 µl
Cell number per Nucleofection® Sample	2.5 x 10 ⁵ cells (Lower or higher cell numbers may influence transfection results)	5 x 10 ⁴ cells (Lower or higher cell numbers may influence transfection results)

Table 3: Contents of one Nucleofection® Sample and recommended program

	100 µl Single Nucleocuvette®	20 µl Nucleocuvette® Strip
	2.5 x 10 ⁵	5 x 10 ⁴
pmaxGFP® Vector	_2 μg	_0.4 μg
plasmid DNA (in H ₂ O or TE)	_ 2 μg	_0.4 – 0.8 μg
siRNA	30 – 300nM siRNA (3 – 30 pmol/sample)	30 — 300nM siRNA (0.6 — 6 pmol/sample)
ıcleofector® X Solution	100 µl	20 µl
	DK-100 (for high efficiency) DJ-100 (for high viability)	DK-100 (for high efficiency) DJ-100 (for high viability)
	plasmid DNA (in H ₂ 0 or TE)	2.5 x 10 ⁵

Table 4: Culture volumes required for post Nucleofection® Steps

	100 μl Single Nucleocuvette®	20 μl Nucleocuvette® Strip*
48-well culture plate	400 µl culture medium l	<u> </u>
96-well culture plate	<u> </u>	110 µl culture medium l
Culture medium to be added to the sample post Nucleofection®	400 μl culture medium l	80 μl culture medium l
* Maximum cuvette volume 200 µl		

Table 5: Recommended volumes for sample transfer into culture plate

	100 μl Single Nucleocuvette®	20 μl Nucleocuvette® Strip*
Culture medium to be added to the sample post Nucleofection®	400 μl culture medium l	80 μl culture medium l
Volume of sample transferred to culture plate	complete sample (use supplied pipettes)	90 µI
* Maximum cuvette volume 200 µl		