

# Non-viral and large-scale delivery of gene editing tools for T-cell manufacturing

Hannah Goldfarb Wittkopf, Frank Morrissey-Wetty, Nadessa Poeppel, Meike Zander, Nicole Kukli, Katharina Stolle, Gina Andretta-Beu, Stefanie Müthel, Valeria Annibaldi and Andrea Toell  
Lonza Cologne GmbH, Nattermannallee 1, 50829, Cologne, Germany

## Introduction

In recent years, non-viral methods for cell engineering have emerged as promising alternatives to viral transduction. Amongst virus-free technologies, electroporation is considered the gold standard, with an increasing presence in clinical trials. Historically, non-viral cell transfection methods, including electroporation, have been reported to suffer with poor scalability and lengthy optimization, sometimes linked to unsatisfactory transfection efficiency at large scale. An ideal technology platform should enable process development and optimization at small scale, with minimal effort required for scaling up and translation to clinical settings.



With our 4D-Nucleofector® LV Unit PRO we provide:

- An optimized cartridge design for improved performance
- Efficient electroporation of complex cargos like CRISPR/Cas9 or Transposon/Transposase
- Easy scale up from 100 µL up to 20 mL
- Electroporation of large volumes with up to  $2.0 \times 10^9$  T cells
- Flexible volumes in fixed-volume cartridge (0.5, 1.0, 1.5 or 2.0 mL)
- Compatibility with cell densities from  $2.5 \times 10^7$  to  $1 \times 10^8$ /mL

## Material and methods

**Cell material:** Cryopreserved human PBMC (Lonza) or CD3+ Pan T Cells (Lonza), activated with TransAct™ (Miltenyi; 3 days) or Dynabeads (Thermo Fisher Scientific®; 2 days)  
**Transfection:** Cells were resuspended in P3 Nucleofector® Solution, cargo was added, the mixture transferred into the required Nucleofector® Vessel and transfected.  
**Analysis:** On day of analysis, TCR alpha (TRAC KO) and/or GFP (KI or transient) expression was evaluated by flow cytometry (NovoCyte, Agilent). Cell count and cell viability were assessed with NucleoCounter® NC-202 (Chemometec).

## Abbreviations

**100 µL** – 100 µL Nucleocuvette® Vessel  
**LV PRO FV** – Fixed volume 2 mL Nucleocuvette® Cartridge PRO  
**LV PRO FT** – Flow-through LV Nucleocuvette® Cartridge PRO  
**LV FV** – Fixed volume 1 mL Nucleocuvette® Cartridge (1st generation)  
**LV** – Large volume

## Transient cargos

### Scalability – mRNA




Figure 1. Transfection efficiency 20 µg/mL eGFP mRNA (n = 2 donors)

### Scalability – plasmid




Figure 2. Transfection efficiency 20 µg/mL plasmid (3.5 kb) or 37.5 µg/mL plasmid (5.0 kb) (n = 2 donors)

## piggybac® Transposon System

### Scalability – activated T cells







Figure 3. Transfection of Transposase mRNA (43 µg/mL) + GenCircle™ Transposon (100 µg/mL; n = 2 donors)

### Scalability – resting T cells



## CRISPR/Cas9 Knock-out

### Scalability – TRAC KO






Figure 4. KO in the constant region of TCRα (TRAC) with 0.77 µM Cas9 and sgRNA in ratio 1:2 (for FV n = 4 donors; for FT n = 1 donor)

### Cell proliferation



## CRISPR/Cas Knock-in

### Scalability – TRAC-GFP KI



### Scalability – Different input volumes




Figure 5. TRAC KI with 0.77 µM Cas9 and sgRNA in ratio 1:2 plus 50 nM dsDNA HDRT (n = 6 donors)

### Scalability – Different HDRT sizes




Figure 6. KI efficiency with 0.77 µM (TRAC) or 0.93 µM (RAB11A) Cas9 and sgRNA in ratio 1:2 plus dsDNA HDRT (n = 2 donors)

### Scalability from 100 µL to 10 mL – TRAC-GFP KI




Figure 7. Editing with 0.77 µM Cas9 and sgRNA in ratio 1:2 plus 50 nM dsDNA HDRT (n = 2 donors)

### Scalability from 100 µL to 20 mL – RAB11A-GFP KI




Figure 8. Editing with 0.93 µM Cas9 and sgRNA in ratio 1:2 plus 8.4 µM dsDNA HDRT (n = 1 donor)

### Robustness – RAB11A-GFP KI




Figure 9. KI efficiency with 0.93 µM Cas9 and sgRNA in ratio 1:2 plus dsDNA HDRT

## Double editing

### TRAC-KO + piggybac® Transposon System




Figure 10. Editing efficiency with Transposase mRNA (86 µg/mL) + GenCircle™ Transposon (200 µg/mL) and TRAC KO (0.77 µM RNP)

### CRISPR/Cas9 Double KO of TRAC and CD45




Figure 11. Efficiency with TRAC KO + CD45 KO for 0.77 µM Cas9 + sgRNA in ratio 1:2 each (n = 2 donors)

## Summary

The next generation 4D-Nucleofector® LV Unit PRO can reliably support non-viral manufacturing of genetically modified T cells. It helps enable:

- Reliable, robust and efficient delivery of complex, clinically relevant cargos
- Easy scale up of the cell engineering process for up to 1 billion cells

Learn more.



Lonza Cologne GmbH – 50829 Cologne, Germany

This research was conducted under a non-exclusive sublicense granted by Arsenal Biosciences, Inc. to U.S. and ex-US patents and pending applications directed to targeted nonviral DNA insertion. Lonza Group Ltd. and its affiliates (collectively and individually, "Lonza") make efforts to include accurate and up-to-date information. However, Lonza makes no representations or warranties, express or implied, including as to accuracy or completeness of information. All trademarks belong to Lonza, and are registered in the USA, EU and/or CH, or used in common law, or belong to third-party owners and are used for only informational purposes. All third-party copyrights have been reproduced with permission from their owners. The user bears the sole responsibility for determining the existence of any third-party rights and obtaining any necessary licenses and approvals. For more information, including regarding legal disclaimers, Lonza's intellectual property rights, and how Lonza collects, uses and protects personal information: [www.lonza.com/legal](http://www.lonza.com/legal), <https://www.lonza.com/about-us/strategy/intellectual-property> and <http://www.lonza.com/privacy>.

©2025 Lonza. All rights reserved.

CT-PO026 09/25