



PyroTec[®] PRO Robotic Solution

FAQ

Move from manual benchtop processes to an automated workflow

Traditional endotoxin testing is plagued by extensive manual steps that increase the risk of human errors, data integrity issues, and inefficient workflows.

The PyroTec[®] PRO Robotic System is an automated, plate-based solution that addresses these pain points and helps your lab with consistent assay preparation and protocol execution, which can improve data integrity and reliability.

This FAQ document aims to provide you with the insights into assay automation and how it can transform a busy manual-based endotoxin laboratory into a streamlined, automated system.

What is PyroTec[®] PRO Robotic System?

PyroTec[®] PRO Robotic System is a fully integrated flexible plate-based solution available to the endotoxin market, designed intrinsically to follow the FDA's Process Analytical Technology and Data Integrity Initiatives. The system combines the speed and reproducibility of a robotic liquid handling platform with Lonza's industry leading WinKQCL[®] Analysis and Detection Software.



How much space is required for the PyroTec® PRO System?

The robot is 145 cm wide (199.39 with PyroWave® extension), 78 cm deep and 87 cm tall. A minimum of 10 cm around the instrument is required for proper function.

Depending on the configuration, the weight of the PyroTec® PRO System ranges from 187 – 200 kg.

Tables are available for purchase from Lonza that will accommodate the PyroTec® PRO System size and weight.

Can you provide an overview of product specifications for the PyroTec® PRO Robotic Solution?

Description	Answer
Throughput per run	42 (21 per plate in duplicate w/ duplicate PPC) with five point standard curve with two negative controls
Sample handling	Yes
Prepares sample dilutions	Yes, 80 (includes sample tubes). All 80 tubes are available to make simple to complex dilution schemes.
Prepares auxiliary dilutions	Yes, for dilutions requiring other than liquid reagent water, i.e. Beta-Glucan blocker
Reconstitution of RSE/CSE	No
Vortex of reconstituted RSE/CSE	No
Prepares endotoxin standards	Yes, flexibility to make up to eight standards per assay (most common is five)
Prepares positive product controls	Yes
Prepares negative controls	Yes
Incubates plate	Yes
Reconstitution of LAL reagents	Yes
Pipettes LAL to plate	Yes
Places plate in incubator	Yes
Barcode scanner	1D barcode; LoadingID on deck automated barcode scanner. 2D barcode; various standard formats supported. Barcodes are entered on the template in WinKQCL® Endotoxin Detection and Analysis Software
Ancillary accessories	Tips, sample tubes, reagent reservoirs, and microplates supplied by Lonza certified down to 0.005 EU/mL
Deck setup time	Fully loaded plate: 5 minutes minimum
Preparation of plate	Fully loaded plate: 40 minutes minimum
Runtime in reader	Approximately 1 hour
Total time to result including preparation and run	1 hour 50 minutes for 21 samples tested neat
Sample test volume	100 µL per well; Minimum sample volume – 600 µL for samples to be diluted; 800 µL for samples to be tested neat
Methods supported	Kinetic chromogenic, kinetic turbidimetric, and Recombinant Factor C

Are there any consumables that need to be purchased to run the PyroTec® PRO Robotic Solution other than standard reagents, dilution tubes, and microplates?

Yes, the following items are available:

00229884	Liquid Handling Filtered Tips
00229888	Trough, Disposable, 100 µL
00309639	Trough, Disposable, 25 mL
00303042	5 mL Sample Tube w/ Cap

Are all the consumables endotoxin-free?

All consumables mentioned above, in addition to 13 x 100 mm tubes (P/N N207), and microplates (P/N 25-340) used in the system, are certified to be endotoxin-free (<0.005 EU/mL).

What is LoadingID?

Analysts can now place samples anywhere in the LoadingID sample runners, and the scanned barcode tells the software where each sample is located on each rack. In addition, the system provides information regarding the detection of duplicates, incorrect samples and missing samples.

The LoadingID unit can be integrated into a new system at time of purchase or by retrofitting an existing PyroTec® PRO System. For retrofitting, please note that the deck layout must be re-optimised and system re-qualification will be required.

What are the power requirements?

- 100-240 Volts
- 50/60 MHz
- 1,200 VA

Which type of liquid handling arm does the PyroTec® PRO System use, liquid or air driven?

The air LiHa technology forms a core element of Lonza's automated solution for endotoxin testing. The air LiHa™ uses air displacement technology to offer increased flexibility and greater convenience for a variety of applications requiring disposable tips. The air LiHa™ has a number of advantages for this application:

- A broad dispensing volume ranging from 10 µL – 800 µL per channel
- A low-level tip ejector that allows removal of disposable tips without risk of cross-contamination due to generation of aerosols
- Utilizes cLLD-capacitive liquid level detection – for fast liquid-level detection and preferential non-contact dispensing, for optimum pipetting performance
- Offers fast dispense down to 10 µL with CVs below 6%
- Each channel is protected by an inline filter

WinKQCL® Endotoxin Detection and Data Analysis Software and Integration

How are manual steps controlled by software?

Pop-up windows in the WinKQCL® Software application guide the user through manual preparation (placement of samples, reagents and consumables on the deck), and every step has to be actively confirmed by the user, which acts as a fail-safe. For data integrity purposes, it reduces the risk of operator error in setting up the deck.

Can manual and automated systems share the same WinKQCL® Software database for trending and data security?

Yes, the system is able to leverage the robust enterprise features of the WinKQCL® Software. This includes storing all raw data and meta data in a WinKQCL® Software database that can be shared amongst the organization.

The software provides fine-grain user access controls and analytics, including trending across the historical records to assess excursions and utilize metadata to assist in root cause analysis investigations.

Can the PyroTec® PRO System interface with LIMS?


Yes, the end user can setup and define the trigger for the data transfer (i.e. import worklists and sample information and export assay data). Lonza offers services to assist with this.

Can plate dilution schemes be transferred to WinKQCL® Software if contained in a LIMS?

Yes, this information can be pulled from a LIMS if available. For example, an export created in CSV can be exported back into LIMS.

Is the PyroTec® PRO System compatible with other, non WinKQCL® Software?

The automation module providing dynamic scripting and control of the robotic system is an integral part of Lonza's WinKQCL® Version 6 Endotoxin Analysis and Detection Software. It is not compatible with any other endotoxin detection software.



Process

Is it possible to automate the standard dilutions?

The standard endotoxin dilutions are created by the system based on what is selected in the WinKQCL® Software Template specific for that assay. The PyroTec® PRO System can dynamically create the dilution scheme in the WinKQCL® Software template.

The system does not reconstitute the stock endotoxin. The endotoxin must be reconstituted and vortexed offline by an analyst prior to running an automated assay – and then placed on the robot deck. The system uses the reconstituted endotoxin to create the standards.

Can you provide more information on the validation of endotoxin standards made using pipetting vs. vortexing?

The only manual step required is for the analyst to perform the initial reconstitution and vortex the initial stock CSE (or RSE) solution prior to running the assay. The robot prepares the standard curve. As we do not place a vortex mixer on the deck, the robot uses multiple aspirate/dispense steps to mix the standards, starting from the stock CSE solution. Study data confirms that this method is comparable to using a standard vortex. This method is included in the reagent package inserts for automation customers.

How does the PyroTec® PRO System mix samples after sample dilution?

All mixing on the system is done by the liquid handling robotic arm described above using similar multiple aspiration/dispense steps to those used in standard preparation. Sample pre-mixing helps ensure heterogenous sample mix prior to dilution.

How does liquid level detection work if different sample containers are used (conductive or pressure based?)

Conductive liquid detection is used by the PyroTec PRO® System. The system currently supports sample handling in 13 x 100 mm tubes (5 mL capped tubes coming soon).

Can you run different samples (different treatment, dilution etc.) on one plate?

The system can accommodate a variety of samples on one plate with different dilutions applied to each sample. The system will apply whatever LAL reagent water dilutions are assigned to each sample on the WinKQCL® Software-generated endotoxin template. Each product

(i.e. sample) defined in the WinKQCL® Software Database is configurable in the automation module software to use a specific liquid class, which determines the pipetting parameters to be applied (different liquid classes for aqueous, viscous, foamy, etc). When a new sample in the WinKQCL® Software Database is detected by the automation software, it is initially assigned a default liquid class, but a different liquid class can be assigned by the user if necessary. The automation software is also able to handle auxiliary dilutions such as beta glucan blocker (P/N N190) and Tris buffer (P/N S50-642). Auxiliary dilutions are configurable per sample. The order in which the auxiliary dilutions are to be applied is also configurable. In addition, the reagent water itself can be configured as an auxiliary dilution if necessary.

Am I limited in the amount and type of sample dilutions?

The user is limited only by the 80 spaces available for sample tubes on the deck layout. The user can create a multitude of sample dilutions with different sample types, utilizing LAL reagent water in addition to auxiliary dilutions simultaneously within the same run.

How many plates can be processed with the instrument?

The system can accommodate up to two independent microplates on a single run of the instrument (i.e. a single deck setup) provided that there are 2 readers associated with the system.

What are the limiting factors (reagents, consumables etc.)?

The limiting factor is deck space on the instrument. Currently a 2-reader instrument system has enough space to accommodate the labware required to run two different types of assays in a single run, supporting chromogenic, turbidimetric, and recombinant factor C assays.

Is it possible to reload the instrument during a run with reagents, consumables, samples? Can you continuously load? Is an automatic retest supported?

After the last microplate is placed in the reader and control is transferred to the WinKQCL® Software to begin reading the plate, the instrument's work is complete. At that time, it would be possible to start preparing the deck for the next run. However, the next run should not be started until all microplate readers (either 1 or 2) have completed reading the microplates, and all microplates have been removed from the readers by the analyst. When preparation is performed on two assays in a single run, the assays are performed sequentially, not concurrently.

The instrument will not start the microplate setup for the second assay until the instrument has placed the microplate from the first assay in the reader to be read by the WinKQCL® Software. When the WinKQCL® Software begins reading the first microplate, the instrument will begin processing the microplate setup for the second assay. The system does not support an automatic retest of a failed sample. The analyst would need to add the sample to the next run to retest it.

How long does the assay take, and can an assay be run overnight?

Typical runtime from start of the automated microplate setup to completion of WinKQCL® Software reading the plate for a single assay is 1.5 hours – 2.5 hours, depending on the number of samples and the complexity of the dilutions. Therefore, a two-assay run will typically complete in 3 hours based on the simplest testing scheme. The system begins setup and preparation of the second assay while the first plate is being measured.

Because the PyroTec® PRO System is automated, once the assay is set up and the run initiated, the run can complete without human interaction, so it may run overnight or between shifts.

Service

How long does it take to perform the installation, operational and performance (IOPQ)?

There are many components that need to be installed and validated. Depending on the customizations needed, estimated time for installation and operational qualification is one week.



Can you provide more information on the IOPQ?

1. Installation Qualification (IQ) – ensures that the hard-ware and software are installed and configured properly in preparation for further testing. Basically, this verifies that the software is installed and that readers and incu-bators have been configured correctly.
2. Operational Qualification (OQ) – ensures that hardware and software components are in proper working order (confirms functionality of LiHa and RoMa arms and functionality to transfer microplates to/from readers).
3. Performance Qualification (PQ) – evaluates the overall functional performance of the system by performing actual automated endotoxin tests using Lonza endotoxin testing reagents. The technician runs an assay to challenge the repeat pipetting accuracy using a Lonza PQ routine test template with samples and dilutions.

What is involved in the annual Preventive Maintenance?

Annual Preventative Maintenance (PM) is included in the Te-Care™ service contract, offered through the Lonza-Tecan partnership. The scheduled PM includes general cleaning, greasing, adjustments and replacement of wear-and-tear parts. Any required updates for instrument reliability, usability and safety will be provided.

Note that Te-Care™ also provides access to Technical Support through Lonza's Scientific Support team, and on-site repairs.

What additional maintenance activities are required?

Routine maintenance tasks are detailed in the User Manual. They include daily inspection and cleaning of the DiTi cones and an inline filter test. Post use, the worktable should be cleaned and DiTi waste emptied. Weekly, a pressure leak test should be performed to ensure tip accuracy.

Who should I contact for technical support?

Your first line of support is always Lonza Scientific Support (numbers are listed at the end of this document). Lonza's team of dedicated professionals will assess your issue and, if required, obtain further support from Lonza Subject Matter Experts depending on your needs.

Our lab would benefit from utilizing the PyroTec® PRO Robotic Solution. What would we need to do prior to scheduling an install?

The best next step would be to contact your sales representative to schedule a meeting to discuss:

- Sample throughput, workflows and types
- Site requirements
- Service and support
- Return on Investment.

Contact us

North America

Customer Service: +1 800 638 8174 (toll free)
order.us@lonza.com
Scientific Support: +1 800 521 0390 (toll free)
scientific.support@lonza.com

Europe

Customer Service: +32 87 321 611
order.europe@lonza.com
Scientific Support: +49 221 99199 400
scientific.support.eu@lonza.com

International

Contact your local Lonza distributor
Customer Service: +1 301 898 7025
scientific.support@lonza.com

Lonza Walkersville, Inc. – Walkersville, MD 21793

Learn more.



All trademarks belong to Lonza, registered in USA, EU or CH or to third party owners and used only for informational purposes. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights. The user bears the sole responsibility for determining the existence of any such third party rights, as well as obtaining any necessary licenses. For more details: www.lonza.com/legal.

©2025 Lonza. All right reserved.

RT-SP015 07/25

lonza.com/endotoxin-automation