



## β-G-Blocker

# Content

Section	Page No.	
1	Intended Use	1
1	Warning	1
1	Explanation of Test	2
1	Principle	3
1	Reagents Supplied and Storage Conditions	3
1	Procedure	3
1	For Use With Kinetic-QCL™ or PYROGENT™-5000 Reagents	3
2	Recommended Procedure	4
2	Interpretation of Results	4
2	References	5

## Important: Read Entire Brochure Before Performing Test

### Intended Use

$\beta$ -G-Blocker is intended as an adjunct for the *In Vitro* end-product endotoxin testing of human and parenteral drugs, biological products, and medical devices using Limulus Amebocyte Lysate (LAL). This product is not indicated as part of a test system for the detection of endotoxin in clinical samples or as an aid in the diagnosis of human disease.  $\beta$ -G-Blocker blocks the reactivity of LAL to  $\beta$ -1,3-Glucans, conferring increased endotoxin specificity to the LAL test. The Pharmacopeial guidelines for LAL testing allows for the neutralization of inhibiting or enhancing substances if the sensitivity of LAL is unaffected<sup>1</sup>.

### Warning

For *In Vitro* use only. Not for the *In Vitro* determination of Endotoxemia in man. Use this reagent for testing with Lonza LAL Reagent only.

## Explanation of Test

The LAL test is an enzyme-mediated assay for the detection of endotoxin<sup>2,3</sup>. This test has largely replaced the rabbit pyrogen test for endotoxin due to the ease of use, cost, sensitivity and specificity of the LAL test. Several factors have been identified that can cause false-positives in the LAL test, i.e. blood products and polynucleotides<sup>4</sup>. Another important source of false-positives involves a class of substances known as  $\beta$ -1,3-Glucans. These substances contain glucose polymers of varying molecular weight linked primarily through  $\beta$ -1,3 glycosidic linkages. If a sufficient amount of  $\beta$ -1,3-Glucans of a particular size class are present, an LAL response may occur that is independent of the endotoxin-mediated response<sup>5,6</sup>.

$\beta$ -G-Blocker is indicated in cases where  $\beta$ -1,3-Glucan contamination is suspected. Examples include preparations derived from yeast and cellulosic materials, including hemodialysis filters. A common dilution/response pattern is seen in the LAL test with samples contaminated with  $\beta$ -1,3-Glucans. This includes a negative response with concentrated samples, a positive response with increasing dilution and eventual negative response at highest dilution. Additionally, with kinetic methods, a synergistic response (i.e. enhancement) is frequently seen with  $\beta$ -Glucan samples spiked with endotoxin.

## Principle

$\beta$ -G-Blocker blocks the glucan-sensitive G pathway in LAL, rendering the test more specific to endotoxin.

## Reagents Supplied and Storage Conditions

### Reagent (B50-700) Yellow-Labeled Vial

$\beta$ -G-Blocker, sterile-filtered liquid, 5 × 5 ml/vial, endotoxin-free.

Store at 2–8°C. Product is stable for four weeks once opened if not contaminated and stored at 2–8°C.

## Procedure

Follow the LAL test procedure in the package insert for either the Kinetic-QCL™ or PYROGENT™-5000 Assays.

## For Use With Kinetic-QCL™ or PYROGENT™-5000 Reagents

If sample dilution is required, first dilute sample with LAL Reagent Water or other suitable diluent. Dilute sample to twice the normal test concentration. For example, if sample is normally diluted 1/10, only dilute to 1/5. After sample dilution has been performed, add one part  $\beta$ -G-Blocker to one part sample. This effectively dilutes the sample 1/2.

## Recommended Procedure

To obtain a 1:1 ratio of product to blocker, add no less than 250 µl of test product dilution to 250 µl of blocker in a depyrogenated glass vessel. Vortex this solution and dispense 100 µl to appropriate wells of a 96-well plate. A control (test product with no blocker, same dilution) is required. A positive product control is required with both the control and the product/blocker at each dilution. It is also recommended to test the blocker diluted 1:1 with LAL Reagent Water and a positive product control.

## Interpretation of Results

Compare the results of the sample tested with and without blocker. Significantly lower values in samples treated with blocker are indicative of glucan contamination.

**Note:** The presence of glucan in a sample typically acts synergistically with any endotoxin to produce an LAL response higher than either one alone (an enhancement effect). For example, a sample with glucan (but no endotoxin) may yield a value of 0.4 EU/ml. The same sample spiked to 0.5 EU/ml endotoxin may yield a value of 1.5 EU/ml.

## References

1. United States Pharmacopeial Convention. General Chapter <85> Bacterial Endotoxins Test. *United States Pharmacopeia* (USP).
2. Levin, J. and F.B. Bang. Clottable protein in *Limulus*: its localization and kinetics of its coagulation by endotoxin. *Thromb. Diath. Haemorrh.* 19:186–197 (1968).
3. Iwanaga, S. The *Limulus* clotting reaction. *Curr. Opinion in Immunol.* 5:74–82 (1993).
4. Elin, R. and S. Wolff. Nonspecificity of the *Limulus* Amebocyte Lysate test: Positive reactions with polynucleotides and proteins. *J. Infect. Diseases.* 128:349–352 (1973).
5. Morita, T., S. Tanaka, T. Nakamura, and S. Iwanaga. A new (1-3)- $\beta$ -D-glucan-mediated coagulation pathway found in *Limulus* amebocytes. *FEBS Lett.* 129:318–321 (1981).
6. Pearson, F. C., J. Bohon, W. Lee, G. Bruszer, M. Sagona, G. Jakubowski, R. Dawe, D. Morrison, and C. Dinarello. Characterization of *Limulus* Amebocyte Lysate-Reactive Material from hollow-fiber dialyzers. *Appl. Environ. Microbiol.* 48:1189–1196 (1984).

[www.lonza.com/pharmabiotech](http://www.lonza.com/pharmabiotech)

Certificate of Analysis: [www.lonza.com/coa](http://www.lonza.com/coa)

## Contact Information

### North America

Customer Service: +1 800 638 8174 (toll free)

Fax: +1 301 845 8338

[order.us@lonza.com](mailto:order.us@lonza.com)

Scientific Support: +1 800 521 0390 (toll free)

[scientific.support@lonza.com](mailto:scientific.support@lonza.com)

### Europe

Customer Service: +32 87 321 611

[order.europe@lonza.com](mailto:order.europe@lonza.com)

Scientific Support: +32 87 321 611

[scientific.support.eu@lonza.com](mailto:scientific.support.eu@lonza.com)

### International

Contact your local Lonza distributor

Customer Service: +1 301 898 7025

[scientific.support@lonza.com](mailto:scientific.support@lonza.com)

### International Offices

Australia	+61 3 9550 0883
Belgium	+32 87 321 611
Brazil	+55 11 2069 8800
France	0800 91 19 81 (toll free)
Germany	0800 182 52 87 (toll free)
India	+91 124 6052941
Japan	+81 3 6264 0660
Luxemburg	+32 87 321 611
Singapore	+65 6521 4379
The Netherlands	0800 022 4525 (toll free)
United Kingdom	0808 234 97 88 (toll free)

---

Lonza Walkersville, Inc. – Walkersville, MD 21793

Unless otherwise noted, all trademarks herein are marks of the Lonza Group Ltd or its affiliates. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights.

©2020 Lonza. All rights reserved.

08300 PN190-6 04/20

RT-MN006

---