

MODA-EM[®] Module v4.1



Paperless quality and manufacturing

The MODA-EM[®] Module is a compliant paperless solution that automates environment monitoring (EM) and quality processes. Easily manage and report on the full spectrum of EM and quality data with timely, accurate monitoring, and benefit from on-demand reporting, trending, and visualization capabilities for in-depth process analysis and ad hoc queries.

Version 4.1 of the MODA-EM[®] Module now offers full support for Windows[®] 11 and SQL Server, and introduces new cell auditing capabilities for enhanced audit trails. Customers will also benefit from new reports, performance enhancements, and enhanced configuration management to help simplify and streamline their pharmaceutical quality operations while maintaining compliance.

MODA-EM[®] Module v4.1: What is new?

Validated Windows[®] 11 and SQL Server compatibility
MODA-EM[®] Module v4.1 functionality, including upgraded

reports, is now fully compatible with Windows[®] 11 for a continuing seamless experience.

Real-time cell auditing for enhanced audit trails

Automatically track data changes in the Sampling, Testing, and Results Entry screens as data is entered or edited to help more comprehensive audit trails.

Powerful new reports provide greater insights

Gain deeper visibility with a range of new reports:

- **Sample Verification Report:** shares information about samples (from scheduled to completed) based on the selected plan and date range
- **Limits and Associated Plans Report:** helps identify plans with default limits that need to be updated
- **Batch Record Packet by Plan Report:** provides a sampling results summary for a specified Product Batch/Lot, Plan, and Test Method

Enhanced performance

Get faster performance with upgrades to FDC cache data queries, schedule panel loading, and Review/Approve queries.

Flexible batch management options

Manage batches manually or through ERP system integration via the import tool, with easy-to-follow set-up instructions.

Easier enterprise device management

Update settings by device type and without XML files, improving ease and efficiency of device management.

Simplified configuration management

Easily manage EM, FDC, and jobs configurations via the MODA® Service Site instead of individual config files.

Figure 1. The new cell auditing function captures data changes within the Sampling, Testing, and Results Entry screens to build comprehensive audit trails.

Figure 2. New date restriction feature reduces data inputting errors.

Figure 3. Access all settings via the MODA® Service Site aiding flexible batch management.

Figure 4. Add jobs easily via the MODA-EM® Module's simplified configuration management.

Ready for enhanced paperless quality and manufacturing?

Contact the Lonza Informatics Team today to discover more about the benefits of the MODA-EM® Module v4.1 and request a demo: <https://bioscience.lonza.com/moda-demo-request-form>

Contact us

North America

Customer Service: +1 800 638 8174 (toll free)
 order.us@lonza.com
 Scientific Support: +1 800 521 0390 (toll free)
 scientific.support@lonza.com

Europe

Customer Service: +32 87 321 611
 order.europe@lonza.com
 Scientific Support: +49 221 99199 400
 scientific.support.eu@lonza.com

International

Contact your local Lonza distributor
 Customer Service: +1 301 898 7025
 scientific.support@lonza.com

Lonza Walkersville, Inc. – Walkersville, MD 21793

Learn more.



All trademarks belong to Lonza, registered in USA, EU or CH or to third party owners and used only for informational purposes. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights. The user bears the sole responsibility for determining the existence of any such third party rights, as well as obtaining any necessary licenses. For more details: www.lonza.com/legal.