

TheraPEAK® ACK Lysing Buffer

Instructions for use

Product use

TheraPEAK® ACK Lysing Buffer is used to lyse red blood cells in a preparation that contains white blood cells.

Formula

NH ₄ Cl	8,024 mg/L
KHCO ₃	1,001 mg/L
EDTA·Na ₂ ·2H ₂ O	3.722 mg/L

Protocol

TheraPEAK® ACK Lysing Buffer was designed for the lysing of red blood cells from mouse lymphocyte preps. The buffer also works for human red blood cells. This product is not to be used for lysing whole blood. Please note that this product is not for *in vitro* diagnostic use.

1. Spin blood, remove supernatant.
2. Add equivalent volume of TheraPEAK® ACK Lysing Buffer to the remaining pellet. The volume should be comparable to that of the pellet.
3. Gently resuspend, swirl tube for 30 – 60 seconds.
4. Fill tube with medium that does not contain serum and centrifuge.
5. Again remove supernatant and retain pellet.
6. This process can be repeated one more time if necessary.
7. Proceed with normal cell protocol.

Specifications

Sterility	pH	Osmolality (mOsm)	Endotoxin (EU/mL)
Negative	7.3-7.6	285-300	≤0.25

Storage

15°C to 30°C

Ordering information

Catalog number	Description	Size
BP10-548E	TheraPEAK® ACK Lysing Buffer	100 mL

Product use statement

This product is not for human or animal *in vivo* or diagnostic use, including use as a diluent or as an excipient. This product is suitable for GMP manufacturing.

The TheraPEAK® ACK Lysing Buffer is covered by patent and/or patent pending rights owned by the Lonza Group Ltd or its affiliates. All trademarks belong to Lonza Group Ltd. and its affiliates (collectively and individually, "Lonza"), and are registered in the USA, EU and/or CH, or used in common law, or belong to third-party owners and are used for only informational purposes. All third-party copyrights have been reproduced with permission from their owners. User assumes all risks of product use and handling. Lonza makes efforts to include accurate and up-to-date information. However, Lonza makes no representations or warranties, express or implied, including as to the accuracy or completeness of information, or the use or handling of these products. The user is responsible for determining if the products supplied by Lonza and the information and recommendations given by Lonza are (i) suitable for intended process or purpose, (ii) are in compliance with environmental, health and safety regulations, and other regulations in the regions and countries where they are purchased, offered for sale, marketed, sold, and used and (iii) will not infringe any third party's intellectual property rights. The user bears the sole responsibility for determining the existence of any such third-party rights, as well as obtaining any necessary licenses and approvals. For more information, including regarding legal disclaimers, Lonza's intellectual property rights, and how Lonza collects, uses and protects personal information: www.lonza.com/legal, <https://www.lonza.com/about-us/strategy/intellectual-property> and www.lonza.com/privacy.

© 2026 Lonza. CT-OP013 09/26. All rights reserved