

Ham's F12 Medium

Product use

Ham's F12 Medium is a nutrient mixture designed to cultivate a wide variety of mammalian and hybridoma cells when used with serum in combination with hormones and transferrin.

Specifications

Sterility	pH	Osmolality (mOsm)
Negative	7.07-7.40	286-305

Storage

2°C to 8°C

Ordering information

Catalog number	Description	Size
12-615F BE12-615F	Ham's F12 Medium with L-glutamine	500 mL

Product use statement

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* diagnostic procedures.