




COVER ALL TUBING ENDS WITH POLY BAG DUST COVERS.

EACH FINISHED CTC-1323 WILL BE
PACKED IN A 4 MIL. POLY BAG.

TUBING TOLERANCES: 84 IN. OR $> \pm 2.00"$
BETWEEN 60 AND 84 IN. $\pm 1.00"$
BETWEEN 12 AND 60 IN. $\pm 0.50"$
12 IN. OR $< \pm 0.25"$

INSTALL CL-01009 CABLE TIES AT EACH CONNECTION.

BILL OF MATERIALS			*NOT SHOWN IN DRAWING
PART#	Q.P.A.	DESCRIPTION	
S4254	1.00	5 LITER BAG WITH FITMENTS	
BC-0005-H05	1.00	5 LITER DIE HANGING STYLE	
PL-01048	-	PLATINUM ULTRAPAK INNER FILM	
PL-01049	-	PLATINUM ULTRAPAK OUTER BARRIER FILM	
TB-01062	3.00	BAG PORT WITH 1/4" BARB	
PK-01039	1.00	0.325" DIA. HANGING BAR 7.5' LONG	
TB-01130	1.50 FT.	1/4" I.D. X 3/8" O.D. PVC TUBING	
TB-01023	1.50 FT.	1/4" I.D. X 3/8" O.D. PCMGs TUBING	
TB-01254	1.50 FT.	1/4" I.D. X 3/8" O.D. C-FLEX NAD TUBING	
CL-01021	3.00	SMALL PINCH CLAMP	
CN-01055	1.00	FEMALE MPC CONNECTOR WITH 1/4" BARB	
CP-01020	1.00	MALE MPC PLUG	
CN-01065	1.00	MALE MPC CONNECTOR WITH 1/4" BARB	
CP-01068	1.00	FEMALE MPC CAP	
CN-01061	1.00	FEMALE LUER WITH 1/4" BARB	
CN-01225	1.00	SWABABLE VALVE	
CP-01055	1.00	MALE COVER CAP	

REV DATE	DESCRIPTION
A 03/29/06	NEW DEVELOPMENT.
B 04/24/07	COMPANY NAME CHANGE TO LONZA, FILM DESCRIPTION CORRECTED TO PLATINUM ULTRAPAK, NO CHANGE TO DESIGN.

Lonza Walkersville, Inc.
8830 Biggs Ford Road
Walkersville, MD 21793
Phone 301.898.7025

CUSTOMER		APPROVALS	DATE
TITLE		DTM	04/24/07
5 LITER BAG		JWP	
PRICE		PROD.	
—		CUSTOMER	
Scale		MARKETING	
—		Engineering	
DRAWING NO.		REV.	
CTC-1323		B	