

Description	CAS #	Chemical Formula	Concentration		Molarity	
			g/L	mg/L	mM	uM
Ammonium Chloride	12125-02-9	NH ₄ Cl	8.024	8.024E+03	150.000	1.500E+05
Potassium Bicarbonate	298-14-6	KHCO ₃	1.001	1.001E+03	10.000	1.000E+04
EDTA (Ethylenediaminetetraacetic Acid Disodium Dihydrate)	6381-92-6	C ₁₀ H ₁₄ N ₂ O ₈ Na ₂ ·2H ₂ O	3.720E-03	3.720	9.994E-03	9.994