

## Penicillin-Streptomycin

Penicillin-Streptomycin is a broad spectrum bacteriostatic and bacteriocidal, with activity against gram negative and gram positive organisms.

### Usage

The suggested working concentration of Penicillin-Streptomycin is 100 units of Potassium Penicillin and 100 µg of Streptomycin Sulfate per 1 mL of culture media.

Cat. no.	Product	Suggested usage volume
17-603 DE17-603	5,000 units/mL Potassium Penicillin 5,000 µg/mL Streptomycin Sulfate	20 mL/L
17-602 DE17-602	10,000 units/mL Potassium Penicillin 10,000 µg/mL Streptomycin Sulfate	10 mL/L
09-757	20,000 units/mL Potassium Penicillin 20,000 µg/mL Streptomycin Sulfate	5 mL/L
17-718	25,000 units/mL Potassium Penicillin 25,000 µg/mL Streptomycin Sulfate	4 mL/L

### Ordering information

#### Penicillin-Streptomycin

Cat. no.	Product	Description
17-603E DE17-603E	Penicillin-Streptomycin 50x Mixture	Contains 5,000 units Potassium Penicillin and 5,000 µg Streptomycin Sulfate per mL in 0.85% saline, 100 mL
17-602E DE17-602E	Penicillin-Streptomycin 100x Mixture	Contains 10,000 units Potassium Penicillin and 10,000 µg Streptomycin Sulfate per mL in 0.85% saline, 100 mL
17-602F DE17-602F	Penicillin-Streptomycin 100x Mixture	Contains 10,000 units Potassium Penicillin and 10,000 µg Streptomycin Sulfate per mL in 0.85% saline, 500 mL
09-757F	Penicillin-Streptomycin 200x Mixture	Contains 20,000 units Potassium Penicillin and 20,000 µg Streptomycin Sulfate per mL in 0.85% saline, 500 mL
17-745E	Penicillin-Streptomycin- Amphotericin B (100X)	10,000 units/ml Potassium Penicillin, 10,000 µg/ml Streptomycin Sulfate, 25 µg/ml Amphotericin B

### Product use statement

#### THESE PRODUCTS ARE FOR RESEARCH USE ONLY.

Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* diagnostic procedures.

All trademarks herein are marks of Lonza Group or its subsidiaries.