

Dextrose-Gelatin-Veronal (DGV)

Product use

DGV is used as a diluent and/or stabilizer in serological testing.

Formula

CaCl ₂	20 mg/L
MgSO ₄ ·7H ₂ O	120 mg/L
NaCl	8,500 mg/L
Gelatin	600 mg/L
Glucose	10,000 mg/L
Sodium Veronal	380 mg/L
Veronal	580 mg/L

For the complete formulation, visit
<https://bioscience.lonza.com>.

Specifications

Routine characterization of DGV includes testing for sterility, pH, osmolality, and endotoxin.

Storage

15°C to 30°C

Ordering information

Catalog number	Description	Size
10-539B	Dextrose-Gelatin-Veronal (DGV)	500 mL

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* diagnostic procedures.

All trademarks herein are marks of Lonza Group or its subsidiaries.