

UltraSALINE A

Product Use

UltraSALINE A is a HEPES buffered saline used in trypsinization protocols of reduced-serum and serum-free systems and preparation of stock solutions.

Formula

Dextrose	720.64 mg/L
KCl	223.68 mg/L
Na ₂ HPO ₄ ·7H ₂ O	268.07 mg/L
Phenol Red	--
HEPES	7150 mg/L
NaCl	7130 mg/L

For the complete formulation, visit www.lonza.com.

Specifications

Routine characterization of UltraSALINE A includes testing for sterility, pH, and endotoxin.

Storage

2°C to 8°C

Ordering Information

Catalog Number	Description	Size
12-747F	UltraSaline A buffered saline without phenol red	500 ml

Product Use Statement

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* procedures.

All trademarks herein are marks of Lonza Group or its subsidiaries.