

RPMI w/ L-Glutamine BE12-702

Description	CAS #	Chemical Formula	Concentration		Molarity	
			g/L	mg/L	mM	uM
D-Glucose anhydrous	50-99-7	C6H12O6	2.000	2.000E+03	11.101	1.110E+04
NaHCO3	144-55-8	NaHCO3	2.000	2.000E+03	23.807	2.381E+04
L-Glutamine	56-85-9	C5H10N2O3	0.300	300.000	2.053	2.053E+03
KCl	7447-40-7	KCl	0.400	400.000	5.366	5.366E+03
NaCl	7647-14-5	NaCl	6.000	6.000E+03	102.669	1.027E+05
Ca(NO3)2.4H2O	13477-34-4	Ca(NO3)2, 4H2O	0.100	100.000	0.423	423.460
Glutathione reduced	70-18-8	C10H17N3O6S	1.000E-03	1.000	3.254E-03	3.254
L-Arginine.HCl	1119-34-2	C6H14N4O2 • HCl	0.242	241.860	1.148	1.148E+03
L-Asparagine.H2O	5794-13-8	NH2COCH2CH(NH2)COOH • H2O	5.682E-02	56.820	0.378	378.472
L-Aspartic Acid	56-84-8	HO2CCH(NH2)CH2CO2H	2.000E-02	20.000	0.150	150.263
Glycine	56-40-6	HO2CCH2NH2	1.000E-02	10.000	0.133	133.209
L-Glutamic Acid	56-86-0	C5H9NO4	2.000E-02	20.000	0.136	135.962
L-Histidine.HCl.H2O	5934-29-2	C6H9N3O2 • HCl • H2O	2.027E-02	20.270	9.671E-02	96.708
L-Hydroxyproline	51-35-4	C5H9NO3	2.000E-02	20.000	0.153	152.520
L-Isoleucine	73-32-5	HO2CCH(NH2)CH(CH3)CH2CH3	5.000E-02	50.000	0.381	381.185
L-Leucine	61-90-5	HO2CCH(NH2)CH2CH(CH3)2	5.000E-02	50.000	0.381	381.185
L-Lysine.HCl	657-27-2	C6H14N2O2 • HCl	4.000E-02	40.000	0.219	218.998
L-Methionine	63-68-3	HO2CCH(NH2)CH2CH2SCH3	1.500E-02	15.000	0.101	100.536
L-Phenylalanine	63-91-2	HO2CCH(NH2)CH2C6H5	1.500E-02	15.000	9.080E-02	90.805
L-Serine	56-45-1	HO2CCH(NH2)CH2OH	3.000E-02	30.000	0.285	285.442
L-Threonine	72-19-5	HO2CCH(NH2)CH(OH)CH3	2.000E-02	20.000	0.168	167.898
L-Tryptophane	73-22-3	C11H12N2O2	5.000E-03	5.000	2.448E-02	24.482
L-Valine	72-18-4	HO2CCH(NH2)CH(CH3)2	2.000E-02	20.000	0.171	170.794
L-Proline	147-85-3	C5H9NO2	2.000E-02	20.000	0.174	173.717
D-Calcium Pantothenate	137-08-6	C18H32CaN2O10	2.500E-04	0.250	5.246E-04	0.525
Folic Acid	59-30-3	C19H19N7O6	1.000E-03	1.000	2.266E-03	2.266
Myo-Inositol	87-89-8	C6H12O6	3.500E-02	35.000	0.194	194.272
Nicotinamide	98-92-0	C6H6N2O	1.000E-03	1.000	8.189E-03	8.189
Pyridoxine HCl	58-56-0	C8H11NO3	1.000E-03	1.000	4.864E-03	4.864
Riboflavine	83-88-5	C17H20N4O6	2.000E-04	0.200	5.314E-04	0.531
Thiamine.HCl	67-03-8	C12H18N4OSCl2	1.000E-03	1.000	2.965E-03	2.965
P-Aminobenzoic Acid	67-03-8	C12H18N4OSCl2	1.000E-03	1.000	2.965E-03	2.965
D-Biotin (Vitamin H)	58-85-5	C10H16N2O3S	2.000E-04	0.200	8.187E-04	0.819
Vitamin B12	68-19-9	C63H88CoN14O14P	5.140E-06	5.140E-03	3.792E-06	0.004
Choline Chloride	67-48-1	HOCH2CH2N(CH3)3Cl	3.000E-03	3.000	2.149E-02	21.485
Phenol Red Salt Sodium	34487-61-1	C19H14O5S	5.310E-03	5.310	1.411E-02	14.107
L-Tyrosine.2Na.2H2O	122666-78-9	C9H9NO3Na2 • 2H2O	2.883E-02	28.830	0.110	110.379
L-Cystine.2HCl	7048-04-6	C3H7NO2S • HCl • H2O	6.519E-02	65.190	0.371	371.241
Na2HPO4 anhydrous	7558-79-4	Na2HPO4	0.800	800.490	5.637	5.637E+03
MgSO4 anhydrous	7487-88-9	MgSO4	4.883E-02	48.830	0.406	405.666