

Grace's Insect Medium

Product use

Grace's Insect Medium is used to cultivate insect cell lines and organs for the propagation of entomopathic viruses.

Product use statement

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* diagnostic procedures.

Specifications

	Sterility	pH	Osmolality (mOsm)	Endotoxin (EU/mL)
04-457	Negative	5.9-6.2	340-380	FIO
04-649	Negative	5.7-6.2	381-411	FIO

Storage

2°C to 8°C

Ordering information

Catalog number	Description	Size
04-457F	Grace's Insect Medium without lactalbumin hydrolysate, yeastolate, gentamicin or FBS	500 mL
04-649F	Grace's Insect Medium with lactalbumin hydrolysate, yeastolate and gentamicin without FBS	500 mL