

Iscove's Modified Dulbecco's Medium (IMDM)

Product use

Iscove's Modified Dulbecco's Medium (IMDM) is a modified DME containing high glucose (4,500 mg/L), sodium pyruvate, additional amino acids, HEPES buffer (in liquid form), selenium and other components. IMDM can support the growth of a wide variety of mammalian cells.

Storage

2°C to 8°C

Ordering information

| Catalog number | Description | Size |
|----------------------|--------------------------------------|--------|
| 12-722F BE12-722F | IMDM with L-glutamine and HEPES | 500 mL |
| 12-726F BE12-726F | IMDM without L-glutamine, with HEPES | 500 mL |

Product Use Statement

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* diagnostic procedures.