



Lucetta™ Connect Software

Find the Lucetta™ Connect Software [SAAF-1003], for data transfer from the Lucetta™ 2 Luminometer to a PC on the enclosed USB stick.

Software Instructions

A guideline how to use the Lucetta™ Connect Software to upload data to a PC can be found on the USB stick and can be downloaded from our webpage via www.lonza.com/lucetta-connect.

Software Installation

Data upload to a PC requires installation of the Lucetta™ Connect Software on the PC.

Content of USB Stick

- Lucetta™ Connect Software Installer
- Lucetta™ Connect Software Instructions
- Lucetta™ 2 Luminometer Manual

Contact Information

North America

Scientific Support:
+1 800 521 0390 (toll free)
scientific.support@lonza.com

Europe

Scientific Support:
+32 87 321 611
scientific.support.eu@lonza.com

©2019 Lonza. All rights reserved.
CD-KI049 02/19

www.lonza.com/lucetta-connect