

## Datasheet For Trypsin/EDTA (10x)

Catalog #: BE02-007E

Material	CAS	MW g/mol	Quantity
NaCl	7647-14-5	58.44	9000 mg/L
EDTA Na <sub>2</sub> ·2H <sub>2</sub> O	6381-92-6	372.24	2000 mg/L
Trypsin (1:250)			1 700 000 U/L

### Contact information

#### North America

Customer Service: +1 800 638 8174 (toll free)  
Fax: +1 301 845 8338  
order.us@lonza.com  
Scientific Support: +1 800 521 0390 (toll free)  
scientific.support@lonza.com

#### Europe

Customer Service: +32 87 321 611  
order.europe@lonza.com  
Scientific Support: +32 87 321 611  
scientific.support.eu@lonza.com

#### International

Contact your local Lonza distributor  
Customer Service: +1 301 898 7025  
scientific.support@lonza.com

Lonza Walkersville, Inc. – Walkersville, MD 21793

For research use only.

All trademarks belong to Lonza or its affiliates or to their respective third party owners. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights.

©2019 Lonza. All rights reserved.  
CD-DS039 05/19

[www.lonza.com](http://www.lonza.com)