



Automation

# Compatibility Matrix

## For WinKQCL® 6 Endotoxin Detection and Analysis Software

The WinKQCL® Endotoxin Detection and Analysis Software offers support for a variety of versions of operating systems and database servers on both 32-bit and 64-bit architectures. Prior to installation, please install all Microsoft® Updates and Service Packs to ensure compatibility.

**No home or starter editions of the operating systems are supported.**

WinKQCL® Software Versions	5.3	6.0	6.0.2	6.1	6.3
<b>Operating System</b>					
Microsoft® Windows® 7	X	X	X	X	X
Microsoft® Windows® 8.x	X	-	-	-	-
Microsoft® Windows® 10	X	X	X	X	X
Microsoft® Windows® 11	-	-	-	-	X
Microsoft® Windows Server® 2008 R2	X	X	X	X	X
Microsoft® Windows Server® 2012 R2	X	X	X	X	X
Microsoft® Windows Server® 2016	X 5.3.4 only	X	X	X	X
Microsoft® Windows Server® 2019	-	*	X	X	X
<b>Database Server</b>					
Microsoft® SQL Server® 2008	X	-	-	-	-
Microsoft® SQL Server® 2008 R2	X	-	-	-	-
Microsoft® SQL Server® 2012	X	X	X	X	X
Microsoft® SQL Server® 2014	X	X	X	X	X
Microsoft® SQL Server® 2016	*	*	X	X	X
Microsoft® SQL Server® 2017	*	*	X	X	X
Microsoft® SQL Server® 2019	*	*	X	X	X

X = Officially Validated

\* = Not validated and will not work by default, but a workaround can be provided to use this version.

- = Not validated, and we do not recommend

PyroTec® PRO System is only compatible with 64 bit Windows® 10 and 11.

## Contact us

### North America

Customer Service: +1 800 638 8174 (toll free)  
order.us@lonza.com  
Scientific Support: +1 800 521 0390 (toll free)  
scientific.support@lonza.com

### Europe

Customer Service: +32 87 321 611  
order.europe@lonza.com  
Scientific Support: +32 87 321 611  
scientific.support.eu@lonza.com

### International

Contact your local Lonza Distributor  
Customer Service: +1 301 898 7025  
Fax: +1 301 845 8291  
scientific.support@lonza.com

Lonza Walkersville, Inc. – Walkersville, MD 21793

All trademarks belong to Lonza, registered in USA, EU or CH or to third party owners and used only for informational purposes. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights. The user bears the sole responsibility for determining the existence of any such third party rights, as well as obtaining any necessary licenses. For more details: [www.lonza.com/legal](http://www.lonza.com/legal).

©2022 Lonza. All rights reserved.

RT-DS031 05/22

[lonza.com/endotoxin-automation](https://lonza.com/endotoxin-automation)