

MODA-ES® Module v4.0

Enhanced Productivity and Usability Through Electronic Batch Records



The Lonza MODA-ES® Module v4.0 is an electronic batch record solution that helps pharmaceutical customers to expedite batch release by minimizing paper-based processes.

The improved module has a wealth of new features to simplify processes, materials, and item management, as well as enhance approvals, data tracking, and external connectivity.

Pharmaceutical companies now get an increased range of tools and functionality to streamline data capture and prevent errors, enabling swifter delivery of medicines to market.

MODA-ES® Module v4.0: What's New?

New Ways to Streamline Manufacturing



Robust weigh and dispense capabilities

Easily and precisely weigh and transfer both active and inactive ingredients, comparing measurements against your bill of materials.



Enhanced materials management

Better manage and maintain material inventory with automatic material level adjustments based on usage, quality status, and expiry. Reduce errors with built-in warnings and restrictions.



Expression builder

Easily create simple and complex formulas with built-in error proofing, formatting, and reusability.



Data promotion tools

Get right first time system data promotion. Create, validate, and move data to a production environment, error-free and without rework to speed up instrument qualification (IQ).




Figure 1. The module now features robust weighing and dispensing capabilities for greater ease and precision

Simpler Use, Easier Integration

The upgraded module also enhances user-friendliness, simplifies approvals, improves tracking, and strengthens external connectivity.

Seamless execution

A modernized HTML 5 user interface makes use easy and intuitive.

Audit trail enhancements

Relevant history tabs contextualize audit and genealogy data, with changes highlighted for at-a-glance updates.

Improved approval configuration

Control who manages approvals, and when. Take your data with you to reduce duplication.

PGP encryption

Safely connect to external systems.

Ready to redefine your manufacturing with enhanced materials management? Contact the Lonza team to request a demo:

bioscience.lonza.com/moda-demo-request-form

Contact Us

North America

Customer Service: + 1 800 638 8174 (toll free)
 order.us@lonza.com
 Scientific Support: + 1 800 521 0390 (toll free)
 moda.support@lonza.com

Europe

Customer Service: + 32 87 321 611
 order.europe@lonza.com
 Scientific Support: + 32 87 321 611
 moda.support@lonza.com

International

Contact your local Lonza Distributor
 Customer Service: + 1 301 898 7025
 Fax: + 1 301 845 8291
 moda.support@lonza.com

Why choose the MODA-ES® Module for manufacturing?

The MODA-ES® Module v4.0 simplifies the management of processes, materials, and items during product manufacturing, so you can expedite product release. Get user-friendly navigation, a modular and configurable design, and easy integration with third party tools and software — all in a cost-effective solution that delivers greater return on investment.

Lonza Walkersville, Inc. – Walkersville, MD 21793
 All trademarks belong to Lonza, registered in USA, EU or CH or to third party owners and used only for informational purposes. All third party copyrights have been reproduced with permission or belong to Lonza. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights. The user bears the sole responsibility for determining the existence of any such third party rights, as well as obtaining any necessary licenses. For more details: www.lonza.com/legal.