

# Agarose and Compatible Techniques

Recovery Method	SeaKem® LE	SeaKem® GTG™	SeaPlaque™	SeaPlaque™ GTG™	NuSieve™ 3:1	NuSieve™ GTG™	MetaPhor™	SeaKem® Gold	SeaPrep™	I.D.NA™
In-gel Reactions				■		■				
β-Agarase			■	■		■				
Phenol/chloroform			■	■		■				
Recovery Columns	■	■	■	■	■	■	■	■		■
Electroelution	■	■	■	■	■	■	■	■		■
Freeze/Squeeze	■	■	■	■	■	■	■	■		■
Blotting										
Southern <1 kb					■					
Southern >1 kb	■									
Northern <1 kb					■					
Northern >1 kb	■									
Specialty Applications										
Viral Plaque Assays			■	■						
PFGE								■		
Cell Culture			■	■					■	
Encapsulation and Embedding of Cells									■	
DNA Identity Testing										■
Comet Assays	■	■								

