

Phosphate Buffered Saline (PBS)

Product use

Phosphate Buffered Saline serves as a wash solution in a wide variety of laboratory applications.

Formula

	17-516	17-517	BE02-017F
KH ₂ PO ₄	144 mg/L	1,440 mg/L	260 mg/L
NaCl	9,000 mg/L	90,000 mg/L	8120 mg/L
Na ₂ HPO ₄	795 mg/L	7,950 mg/L	933 mg/L

***THESE PRODUCTS ARE FOR RESEARCH USE ONLY.** Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* diagnostic procedures.

Storage

15°C - 30°C

Ordering information

Catalog number	Description	Size
17-516F*	PBS (1X) without calcium or magnesium	500 mL
BE17-516F*		
17-516Q*	PBS (1X) without calcium or magnesium	1 L
BE17-516Q*		
BEBP17-516Q		
17-517Q*	PBS (10X) without calcium or magnesium	1 L
BE17-517Q*		
BE02-017F*	PBS (1X) with calcium and magnesium, with phenol red	500 mL

Product use statements

GMP PRODUCTS ARE INTENDED FOR RESEARCH OR FOR FURTHER MANUFACTURING USE ONLY. This product is not intended for direct therapeutic use in humans