

MEM Eagle Vitamin Mixture (100X)

Product Use

MEM Eagle Vitamin Mixture is used to support cell growth in tissue culture media.

Formula

D-Ca pantothenate	100 mg/L
Choline c hloride	100 mg/L
Folic Acid	100 mg/L
i-Inositol	200 mg/L
Nicotinamide	100 mg/L
Pyridoxine-HCl	100 mg/L
Riboflavin	10 mg/L
Thiamine-HCl	100 mg/L

Usage

10 ml per liter

Specifications

Sterility	pH	Osmolality (mOsm)
negative	4.53-7.05	1-18

Cell growth promotion	Endotoxin (EU/ml)
≥75% of control	FIO

Storage

-10°C to -20°C

Product Use Statement

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* procedures.

Ordering Information

Catalog number	Description	Size
13-607C	MEM Eagle Vitamin Mixture (100X)	50 ml