

Penicillin-Streptomycin-Amphotericin B (100X)

Penicillin-Streptomycin-Amphotericin B is a broad spectrum bacteriostatic and bactericidal agent, with activity against gram negative and gram-positive organisms including yeast and mold.

Formula

Amphotericin B	25 µg/mL
NaCl	8,500 mg/L
Potassium Penicillin	10,000 units/mL
Streptomycin Sulfate	10,000 µg/mL

Usage

10 mL per liter

Specifications

Sterility	pH	Osmolality (mOsm)	Endotoxin (EU/mL)
Negative	3.3-7.2	290-362	FIO

Storage

-15°C to -40°C

Ordering information

Catalog number	Description	Size
17-745E	10,000 units/ml Potassium Penicillin and 10,000 µg/mL Streptomycin Sulfate and 25 µg/ml Amphotericin B	100 mL

Product use statement

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. Not approved for human or veterinary use,

for application to humans or animals, or for use in clinical or *in vitro* diagnostic procedures.

All trademarks herein are marks of Lonza Group or its subsidiaries