



Automation

Compatibility Matrix

For WinKQCL® 6 Endotoxin Detection and Analysis Software

The WinKQCL® Endotoxin Detection and Analysis Software offers support for a variety of versions of operating systems and database servers on both 32-bit and 64-bit architectures. Prior to installation, please install all Microsoft® Updates and Service Packs to ensure compatibility.

WinKQCL® Software Versions	6.0	6.0.2	6.1	6.3	6.4	6.5
Operating System						
Microsoft® Windows® 10	X	X	X	X	X	X
Microsoft® Windows® 11	-	-	-	-	X	X
Microsoft® Windows Server® 2016	X	X	X	X	X	X
Microsoft® Windows Server® 2019	-	-	X	X	X	X
Microsoft® Windows Server® 2022	-	-	-	-	X	X
Database Server						
Microsoft® SQL Server® 2016	*	X	X	X	X	X
Microsoft® SQL Server® 2017	*	X	X	X	X	X
Microsoft® SQL Server® 2019	-	-	X	X	X	X
Microsoft® SQL Server® 2022	-	-	-	-	X	X

X = Officially Validated

* = Not validated and will not work by default, but a workaround can be provided to use this version.

- = Not validated, and we do not recommend

PyroTec® PRO System is only compatible with 64 bit Windows® 10 and 11.

No home or starter editions of the operating systems are supported.

Contact us

North America

Customer Service: +1 800 638 8174 (toll free)

order.us@lonza.com

Scientific Support: +1 800 521 0390 (toll free)

scientific.support@lonza.com

Europe

Customer Service: +32 87 321 611

order.europe@lonza.com

Scientific Support: +49 221 99199 400

scientific.support.eu@lonza.com

International

Contact your local Lonza Distributor

Customer Service: +1 301 898 7025

Fax: +1 301 845 8291

scientific.support@lonza.com

Lonza Walkersville, Inc. – Walkersville, MD 21793

Lonza Group Ltd. and its affiliates (collectively and individually, “Lonza”) make efforts to include accurate and up-to-date information. However, Lonza makes no representations or warranties, express or implied, including as to accuracy or completeness of information. All trademarks belong to Lonza, and are registered in the USA, EU and/or CH, or used in common law, or belong to third-party owners and are used for only informational purposes. All third-party copyrights have been reproduced with permission from their owners. The user bears the sole responsibility for determining the existence of any third-party rights and obtaining any necessary licenses and approvals. For more information, including regarding legal disclaimers, Lonza’s intellectual property rights, and how Lonza collects, uses and protects personal information: www.lonza.com/legal, <https://www.lonza.com/about-us/strategy/intellectual-property> and www.lonza.com/privacy.

©2026 Lonza. All rights reserved.

RT-DS031 04/26

lonza.com/endotoxin-automation